

Ergonomically shaped High backrest - Available on our Active Seating System

Easy Reach Braking System, ensuring user safety & stability when getting in and out of the chair



VELA Tango 510

Ability to easily manoeuvre the chair around the office/home

Large tilt backwards :: Stable fold-back armrests :: possibility for fold-back footrest :: Telescopic heightadjustment of backrest









VELA :: Gøteborgvej 8-12 :: 9200 Aalborg SV :: Denmark :: Tel: +45 96 34 76 00 :: mail@vela.dk

VELA Bariatric chair

The VELA Bariatric chair provides an ergonomic sitting posture, the possibility to move around in the home and the ability to be safe and stable due to brakes and armrests

The VELA Bariatric chair is possible to adapt, eg. seat and backrest to fit your needs REINFORCED ARMRESTS ADJUSTABLE ARMRESTS STABLE SEAT WITH POSSIBILITY FOR CUSTOMISATIONS A CENTRAL BRAKE WHICH PROVIDES SAFETY WHEN. GETTING INTO AND OUT OF THE CHAIR ELECTRIC HI/LO LIFT **FUNCTION HEIGHT** ADJUSTMENT, PROMPTING USER INDEPENDENCE STABLE 4 WHEELED BASE · ·

VELA Tango 300El

Ability to manoeuvre the chair around in the home

Secure and reinforced, stable armrests :: Reinforced backrest :: Reinforced base with wheels and brake :: Electric height adjustment









VELA Mobile perching stool

Ergonomically shaped seat and backrest - wheels with brake for optimum use & flexibility

Easy Reach Braking System, ensuring user safety & stability when getting in and out of the chair







VELA :: Gøteborgvej 8-12 :: 9200 Aalborg SV :: Denmark :: Tel: +45 96 34 76 00 :: mail@vela.dk



VELA chairs are developed so that individuals with disabilities can maintain independence and stay active. The right accessories on the chairs can be crucial for providing, e.g. good sitting posture and improved participation for the user.








BODY SUPPORT

- :: Lateral supports are suitable for users who require stability at the sides and front
- :: Thigh/side supports have padded cushions that offer stability and comfort without limiting activity
- :: Abduction blocks help to keep the user in an optimal comfortable position










FOOT- & LEG REST

- :: Footrest with fold-out LP plate. This footrest is fitted to the seat, and travels with the seat when it is height-adjusted
- :: Fold-out foot bar which can be folded up so that it is not in the way or can be used as a step when the user gets into and out of the chair
- :: A foot and ankle harness is suitable for people who find it difficult to keep their feet on the footplate



BELTS

- :: The Y seat belt is a 4-point lap belt with adjustable length. The seat belt is available with either a press lock or a binding lock
- :: The H seat belt controls the users body and provides stability in the sitting position
- :: The Pelvic four-point lap belt is suitable for individuals who require a very supportive lap belt



MISCELLANEOUS

- :: The head support with pad is adjustable in height, depth and angle so the user can achieve optimum support for the head and neck
- :: The multi-adjustable push bar takes account of the healthcare worker's work posture when having to transport a user in the chair
- :: The table tray folds down and can be removed. It is height, depth and width adjusted to match the user's sitting position
- :: Different armrest pillows for optimum support and for providing an ergonomic sitting posture
- :: An extension grip is an advantage when having difficulties reaching a standard length brake stick. A T-brake lever or a V-grip makes it easier for people to use the chair.

ES backrest with lateral O seat/backrest when E seat with pommel support using special cushion

UPHOLSTERY AND FOAM

- :: VELA offer individual solutions within padding, comfort and support on all seats and backrest, e.g. a pommel in front of seat or lateral support
- :: If a specially customised cushion is required, an O seat/backrest can be used as a base to put the customised cushion onto
- :: VELA also solve special needs within Memory foam, Dunlopillo or other solutions



Promoting independence at school & home.

VELA Blues - Powerchairs that 'grows' with the child



THE MODULAR
DESIGN ENABLES
THE SEAT PAD
& BACKREST TO
BE CHANGED TO
ACCOMMODATE
GROWTH



"I like my chair because I can move around with my friends and chase my sister and brother." Sadie Denman "Sadie's posture in the Vela chair is suitable for eating, drin-king, using her communication aid and of course, driving the chair! This makes her the best she can be at school and gives her some much needed independence." Nicola Denman

- :: Great for playground and classroom activities
- :: Compact and manoeuvrable
- :: Easy to get into and out of
- :: Ergonomically-shaped seat and back

Typical users:- MS, MD, CP, fibromyalgia, spinal surgery, chronic musculoskeletal pain etc Finance Options Available

Contact us for further information or to arrange free product training/CPD session Free Assessment Available on Request

VELA :: Gøteborgvej 8-12 :: 9200 Aalborg SV :: Denmark :: Tel: +45 96 34 76 00 :: mail@vela.dk www.vela.eu



Promoting independence at at Work VELA Blues 100 - Compact indoor powerchair, with electric hi-lo

COMPACT
POWERCHAIR
THAT IS EASY
TO MANOUVRE
AROUND THE
OFFICE

marving live





VELA Blues 100

VELA Blues 100-SW



"I have regained my independence and self respect!."

Hayley Rush

"Since receiving my VELA Blues Powerchair, my work life has been made that much easier. I can get myself to the bathroom, the drinks machine and other parts of the office. I have regained my independence and self respect! I no longer have to wait for people to push me around – it's fantastic!! " Hayley Rush

:: Great for Office or Home based Activities

:: Compact and manoeuvrable

:: Easy to get into and out of

:: Ergonomically-shaped seat and back


Typical users:- MS, MD, CP, fibromyalgia, spinal surgery, chronic musculoskeletal pain etc Finance Options Available

Contact us for further information or to arrange free product training/CPD session Free Assessment Available on Request

VELA :: Gøteborgvej 8-12 :: 9200 Aalborg SV :: Denmark :: Tel: +45 96 34 76 00 :: mail@vela.dk



VELA chairs are developed so that individuals with disabilities can maintain independence and stay active. The right accessories on the chairs can be crucial for providing, e.g. good sitting posture and improved participation for the user.









Lateral support

Thigh/side support

BODY SUPPORT

- :: Lateral supports are suitable for users who require stability at the sides and front
- :: Thigh/side supports have padded cushions that offer stability and comfort without limiting activity
- :: Abduction blocks help to keep the user in an optimal comfortable position







Snap Hook



Foot and ankle harnesses

FOOT- & LEG REST

- :: Footrest with fold-out LP plate. This footrest is fitted to the seat, and travels with the seat when it is height-adjusted
- :: Snap Hook for bags is suitable for shopping bags and school bags etc. The snap hook for bags is fitted to the rear of the chair's back frame.
- :: A foot and ankle harness is suitable for people who find it difficult to keep their feet on the footplate



Y seat belt



H seat belt



Pelvic four-point lap belt

BELTS

- :: The Y seat belt is a 4-point lap belt with adjustable length. The seat belt is available with either a press lock or a binding lock
- :: The H seat belt controls the users body and provides stability in the sitting position
- :: The Pelvic four-point lap belt is suitable for individuals who require a very supportive lap belt



Head support



Covers



Table tray



P&GVR2



P&GR-net

MISCELLANEOUS

- :: The head support with pad is adjustable in height, depth and angle so the user can achieve optimum support for the head and neck
- :: Fabric top covers for spillage or water-proof covers for incontinence. A cover is hygienic and can easily be removed and washed
- :: The table tray folds down and can be removed. It is height, depth and width adjusted to match the user's sitting position
- :: P&G VR2 Control Unit is a simple, user-friendly joystick with the option of programming all drive parameters. Available with and without controls for lights and indicators.
- :: P&G R-net unit provides for optimal programming of the electric wheelchair. All drive parameters can be adapted to meet the user's needs. Advanced control and installation of the Bluetooth module for operation of PCs and smartphones etc.



E seat with pommel



ES backrest with lateral support



O seat/backrest when using special cushion

UPHOLSTERY AND FOAM

- :: VELA offer individual solutions within padding, comfort and support on all seats and backrest, e.g. a pommel in front of seat or lateral support
- :: If a specially customised cushion is required, an O seat/backrest can be used as a base to put the customised cushion onto
- :: VELA also solve special needs within Memory foam, Dunlopillo or other solutions