



## More than 30 years at eye level

R82 produces technical aids and appliances for disabled children and young people throughout the world. Based in Denmark with sales companies in England, Holland, Italy, Germany, Australia and the USA, we develop products tailored to individual needs.

Through comprehensive development R82 is constantly moving foward- but without compromising quality of service and closeness to the customer. Distributors in 45 countries spread over most of the globe ensure forward personal support, even after the product has been delivered.

The objective of each and every R82 employee is: adjustment to individual needs to provide the child the greatest possible freedom – in a world of movement.

## Š.

### **SEATING**

<b>1</b>	Kudu7
<i>,</i>	Cougar
	x:panda
	Panda Futura25
	Panda Futura 5
	Scallop
	Wombat Living
	Cricket
	Stingray
	Cruiser
	EZ Rider65
	Rodeo71
	High-low Frame77
	High-low:x Frame81
	Multi Frame
	Multi Frame:x93
	Stingray Frame99
	Combi Frame:x
	Strong Base
	STANDING
	Meerkat
''	Rabbit Up121
	Toucan
	Caribou
	Gazelle PS
'	WALKING
į į	Mustang
	Crocodile
	Pony
	TRANSPORTATION
<u>எ</u>	Stingray Swing-out
	Juligray Jwilig-Out
	TOILET & BATH
	Flamingo High-low
39	Flamingo
	Heron
	Swan
	Manatee
	Orca Penguin





Head Supports	9
Side Supports	3
Vests and belts	7
Travs	1



Hi and welcome to the R82 catalogue for 2019/2020. I hope you like it...





## Seating, overview

(Wheelchairs. wo	systems orking chairs, floor seats)			Kudu	Cougar	Scallop	Wombat Living
Target users	Max. load			60 kg	100 kg	70 kg	70 kg
	Max. load for transportation			57 kg	75 kg	-	-
	User height			850-1700 mm	1300-1700 mm	600-1400 mm	750-1900 mm
Features	Self centering seat plane angle adjust	ment/Back recline		•	•		•
	Safe transportation in vehicles			•	•		
	Back recline			•	•		•
	Seat plane angle adjustment			•	•		•
	Easy Fold - Light weight					•	
Environment	Indoor			•	•	•	•
	Outdoor			•	•	•	
Seating s	systems		Cricket	Stingray	Cruiser	EZ Rider	Rodeo
Target users	Max. load		50 kg	40 kg	113 kg	77 kg	63 kg
	Max. load for transportation		43 kg	40 kg	113 kg	77 kg	63 kg
	User height		850-1400 mm	750-1500 mm	600-1800 mm	1000-1800 mm	700–1700 mn
Features	Safe transportation in vehicles		•	•	•	•	•
	Back recline		•	•			•
	Seat plane angle adjustment			•			•
	Easy Fold - Light weight		•		•	•	•
	180° turnable seat			•			
 Environment	Indoor		•	•	•	•	•
	Outdoor		•	•	•	•	•
Seats					x:panda	Panda Futura	Panda Futura 5
Target users	Max. load				80 kg	60 kg	85 kg
	Max. load for transportation				75 kg	57 kg	85 kg
	User height				650-1800 mm	1000-1400 mm	1600 mm
Features	Dynamic back				•		
	Growth (width, depth, back height)				•		
	Back recline (backrest angle)				-5°-30°	•*	-5°-30°
	Fits on a variety of frames				•	•	•
Environment	Indoor				•	•	•
	Outdoor				•	•	•
		High-low & High-low:x	Multi Frame	Multi Frame:x	Stingray Frame	Combi Frame:x	Strong Base
Frames			4001.7	100 kg	40 kg	100 kg	120 kg
	Max. load	60-100 kg	100 kg				
	Max. load Max. load for transportation	60-100 kg 60 kg*	100 kg	100 kg	40 kg	100 kg	-
Target users Features		_		100 kg	40 kg	100 kg	•
Target users	Max. load for transportation  Height adjustment (manual, gas,	60 kg*		100 kg	40 kg	100 kg	•
Target users	Max. load for transportation  Height adjustment (manual, gas, hydraulic, spindle or power)	60 kg*	100 kg				•
Target users	Max. load for transportation  Height adjustment (manual, gas, hydraulic, spindle or power)  Seat plane angle adjustment	60 kg*	100 kg				•
Target users Features	Max. load for transportation  Height adjustment (manual, gas, hydraulic, spindle or power)  Seat plane angle adjustment  180° turnable seat unit	60 kg*	100 kg				•
Target users Features	Max. load for transportation  Height adjustment (manual, gas, hydraulic, spindle or power)  Seat plane angle adjustment  180° turnable seat unit Indoor	60 kg*	100 kg				• • • 600 mm

4 R82 | a part of Etac



# **b.** Seating

Comfortable seating systems giving the greatest possible support and mobility!



## **Kudu**<sup>®</sup>

The Kudu is a practical wheelchair which combines form and function to provide the right fit and support for a comfortable posture.

The Kudu is available in 4 sizes.

## **Highlights about Kudu**



### **Growth adaptable**

The base of the seat and the back rest are equipped with a mechanism that makes it possible to adjust height and depth as well as mount fix-locks for vests and belts.



Accessories are easily added and removed as needed for maximum user comfort



## **Centre of gravity**

To adapt to different needs, the seat has a centre of gravity making it possible to adjust the rear wheels in 3 different positions.



### Safe transportation

The Kudu complies with ISO 7171-19 for transporting children up to 57 kg in motor vehicles.

Read more at R82.co.uk Seating | Kudu





### **Standard**

Swing-away anti tip device, tilt in space, back angling (recline).

## **Options**

### Size

- Size 1
- Size 2
- Size 3
- Size 4

### Cushions

- Kudu cushions, upcharge
- No cushions
- Cushions with contour, upcharge

### For the hands

- Push brace, angle adjustable
- Push brace, height and angle adjustable
- Push handles, height adjustables

### Foot support

- Foot plate, centrally mounted
- Parted foot rests (available for size 2+3+4)

### Rear wheels

- 12½"/16" solid or air
- 20"/22"/24" solid or air with quick release

### Front wheels

- 6" solid front wheels
- 7" solid or air front wheels

#### **Brakes**

- Wheel brakes (available for 20", 22" and 24")
- Drum brakes

## Anti-tips & spoke protectors

- Anti-tips
- Spoke protectors

### Frame colour







White





Red



Kudu with Push brace, angle adjustable



Kudu with Push handles, height adjustable



## **Measurements**

		1	2	3	4
	Effective seat width	260 mm	300 mm	340 mm	380 mm
<b>₽</b> °	Effective seat depth	220-320 mm	270-370 mm	310-410 mm	350-450 mm
<u>1</u> 2°	Backrest height	370-460 mm	400-520 mm	430-580 mm	430-630 mm
<u>921</u>	Seat surface height at front edge	480 mm	480 mm	480 mm	480 mm
8	Max. load	40 kg	50 kg	60 kg	60 kg
8	Max. load at transportation	40 kg	50 kg	57 kg	57 kg
P	Backrest angle	-5°- 30°	-5°- 30°	-5°- 30°	-5°- 30°
B.	Seat plane angle	-3°- 45°	-3°- 45°	-3°- 45°	-3°- 45°
<u>p</u> ?	Length	690 mm	690 mm	740 mm	740 mm
<u> </u>	Overall width	590 mm	630 mm	670 mm	710 mm
<u>G</u> ?	Folded length	825 mm	825 mm	875 mm	875 mm
<u> = </u>	Folded width	610 mm	650 mm	690 mm	730 mm
ig,	Folded height	625 mm	625 mm	625 mm	625 mm
<b>∂</b> 8▲	Total mass	19 kg	21 kg	23 kg	25 kg
				Further m	easurements at our website

The Kudu is measured with the following wheel combination, placed in the standard wheel position: 22" rear wheel and 6" front wheel.

## Nice to know

The many different configuration possibilities and accessories to the kudu makes it ideal to situations where complex and individualized seating solutions is needed and reusing of the product is desirable. Examples of individualized solutions is placing the centrally mounted foot plate underneath the seat or using the angle adjustable armrest and tray allowing the tray to be horizontally when the seat plate is angled.

The focus of the seating solution will often be comfort, however when configured with large driving wheel and cambering it is possible to promote self-propelling of the wheelchair over short distances to support user independence.

Read more at R82.co.uk Seating | Kudu





## **Accessories**



**Head support** 99xxxxx-xx



**Back extension cushion** 95743-xxx



**Fixed side supports** 95401-xx



Swing-away side supports 8310xx-xxx



Back wedges E6301x-60



**Arm rests/Tray** 95712-x/9576351 **Padding** 95712-x



**Hip supports** 95713-x0



Pommel 95430-x0 Padding 95430-x0260



**Vest** 89113-x



Mud guards 9520xx



**Transport fittings** 95900-10



Adductors 95429-10

For more accessories - visit our website









Pommel 95430-x



**Upper arm supports** 8910407-x0



Angle adjustable arm rests, incl. side plates 95712-x0





**Canopy** 87401x0 (Requires 9589100)





Read more at R82.co.uk Seating | Kudu



## Cougar®

The Cougar is built for rehabilitation with comfort provided as standard and it has many adjustment possibilities and wide range of accessories.

The Cougar is available in 4 sizes for use by children and teenagers.

## **Highlights about Cougar**



### **Back recline**

The adjustment of the back tilt occurs in the anatomical pivot point of the hip and therefore with minimal shear.



### **Tilt in space**

The Cougar has self-centering centre of gravity making it a stable comfortable chair. The tilt in space function gives the option of changing the users position and reduces pressure in the seating area.



### **Height adjustable**

The seat to floor height is adjustable. Adjust the height on the back of the chair and the angle of the front wheels will remain in the correct position.



13

### **Accessories**

Easily interchangeable for user comfort.

Read more at R82.co.uk Seating | Cougar





### **Standard**

Patented change of seat to floor height, prepared for other manufactures seats and backs, easy change of seat depth, selfcentring centre of gravity, patented recline to ensures a minimum of shear, anti-tip safety.

### **Cushion colours**



Black

### Frame colour



Silver

## **Options**

### Four sizes

- Size 1
- Size 2
- Size 3
- Size 4

### Rear wheels

- 12½" solid or air
- 20" solid or air
- 22" solid or air
- 24" solid or air

#### **Brakes**

- Wheel brake (not available for 12½")
- Drum brake

#### Front wheels

- 6" solid or air
- 8" solid or air (not in combination with 20" rear wheels)

### Leg rests

- No leg rests
- Parted leg rests (length adjustable) with gas spring
- Parted leg rests, fixed

### **Calf** support

- No calf support
- With calf support

R82 | a part of Etac

### Seat

- Wooden seat, no cushion
- Seat cushion, contoured
- No seat cushion and no wooden seat board

### Back rest

- Strap back without cushion
- Strap back with thin cushion
- Strap back with contoured cushion
- Wooden back without cushion
- Wooden back with thin cushion
- Wooden back with contoured cushion

#### Back recline

- No recline
- Recline

#### Recline

- Operated from push brace or on side back
- User operated back adjustment
- User operated



Standard Cougar mounted with drum brake and wheel brakes.



Standard Cougar mounted with drum brake, wheel brakes and cushions.



## **Measurements**

		1	2	3	4
	Effective seat width	350 mm	400 mm	450 mm	500 mm
<b>₽</b> °	Effective seat depth	350 - 450 mm	370 – 470 mm	430 – 530 mm	490 – 590 mm
<u>P</u>	Backrest height	450 mm	500 mm	600 mm	600 mm
<u>921</u>	Seat surface height at front edge	460 – 530 mm	460 – 530 mm	460 – 530 mm	460 - 530 mm
8	Max. load	75 kg	75 kg	100 kg	120 kg
8	Max. load at transportation	75 kg	75 kg	100 kg	120 kg
Di	Backrest angle	37°	37°	37°	37°
$\mathscr{Z}'$	Seat plane angle	90°	90°	90°	90°
<u>P</u> ?	Length	780 mm	780 mm	830 mm	830 mm
<u> </u>	Overall width	570 mm	620 mm	670 mm	720 mm
<u>G</u> ,	Folded length	780 mm	780 mm	830 mm	830 mm
<u> - </u>	Folded width	570 mm	620 mm	670 mm	720 mm
ie,	Folded height	960 mm	990 mm	1080 mm	1100 mm
<u>B</u> 8₩	Total mass	25 kg	26 kg	27 kg	28 kg
				Further me	easurements at our websit

The Cougar is measured with cushion and the following wheel combination, placed in the standard wheel position: 22" rear wheel and 6" front wheel. The Cougar is measured with the seat in horizontal and back in vertical position. We supply various types of leg rests and cushions, which is why the weight is listed without. Different settings may result in another set of measurements.

15

## Nice to know

Cougar can be fitted with a wooden board back or a strap back. Both backs are prepared for velcro fit cushion systems. Furthermore, the Cougar is often used for other seats and backs, such as moulded seats.

The stepless adjustment of the vertical position of the side supports allows a precise position of the supports in relation to the need of the individual user.

Read more at R82.co.uk Seating | Cougar





## **Accessories**





**Head supports** 99xxxxx-xx



Back extension, fixed 970055-x13



**Back extension** 970050-x13



Pommel 95430-x



Swing away sidesupport 83xxx-xxx



Arm rest/Extension 970600-x/970630-1



Tray 970240-x



89661-4P



Hip support 970751-x



Leg rest fixed with calf support 970650-x



Spoke protector 97092x



Mudguards 970560

For more accessories - visit our website





Transport fittings 970900




Vent tray 970815



Arm rest with side plates 970600-x Without side plates 970601-x



Leg rest with gas strut



Bar for foot support 897526x Bar for parted foot support 897528x (not shown)



**Recline add-on** 97094x



**Foot plate model Wombat** 95145-2, 95145-3, 95145-4



17

**Parted foot supports** 89146-x



## Freedom of configuration

The Cougar comes with different types of leg rests and foot plates, to suit any user's needs.

Read more at R82.co.uk Seating | Cougar



## x:panda®

x:panda is a multi-adjustable seating system that can be configured to meet a range of postural requirements and to accommodate growth. Its dynamic back allows a child to extend and push back against relatively light resistance. It is available in 4 sizes.

## Highlights about x:panda



The x:panda has a unique dynamic back that absorbs the energy of extension



### Multi-adjustable

Overlap in sizing means that the starting option always offers the maximum growth potential coupled with individual leg positioning to offer an improved base of support and additional support for maintaining the correct degree of abduction/adduction if required.



### Growth

Individual growth in width, depth and back height extends the life of the seat and accomodates the greatest range of sizes for any chair on the market.



### **Back recline**

The patented articulation of the back rest has been designed to eliminate shear and sliding during extension. This means that support remains uniquely in position during movement improving functionality and comfort for the user.

19

Read more at R82.co.uk Seating | x:panda





## **Standard**

Hip supports - front and back with cover, ad-/abduction support with cover, spacers for seat extension, dynamic backrest, lumbar support (x:panda size 4 only), fix lock fittings for mounting of various belts and vests, lockable dynamic gas spring.

## **Options**

### Four sizes

- Size 1
- Size 2
- Size 3 - Size 4

### Seat width

- Small - Madium
- Large

### Back cushion

- Small
- Medium - Large

### Seat cushion

- Small
- Medium - Large

## Seat cushion without abduction

- Small
- Medium
- Large

### **Gas spring**

- Non-dynamic gas spring
- Lockable dynamic gas spring (standard)

It is possible to combine the different seat and back sizes: Size 1 seat with size 2 back, Size 2 seat with size 3 back, Size 3 seat with size 4 back, Size 2 seat with size 2 back, Size 4 seat with size 3 back.

### **Cushion colour back and seat**



Black



Pink\*









Blue Lime green

\*Not available for size 4

### **Cushion colour seat without abduction**



Black







x:panda with abduction support and red cushions

20 Seating | x:panda R82 | a part of Etac



## **Measurements**

		1	2	3	4
	Effective seat width	160, 200, 240 mm	200, 240, 280 mm	260, 310, 360 mm	340, 400, 460 mm
<u>_</u>	Effective seat depth	160-260 mm	220-330 mm	260-390 mm	380-530 mm
Ī	Backrest height	270-370 mm	320-420 mm	350-470 mm	500-710 mm
	Max. load	30 kg	30 kg	60 kg	80 kg
	Max. load at transportation	30 kg	30 kg	60 kg	75 kg
9	Backrest angle	-5°- 30°	-5°- 30°	-5°- 30°	-5°- 30°
	Overall width	310, 350, 390 mm	380, 420, 460 mm	470, 520, 570 mm	590, 650, 710 mm
	Total mass	6,5 kg	8,5 kg	11 kg	13 kg
				Further me	easurements at our website

## Nice to know

The back height, seat width and depth adjustments, different accessories, the 16 pre-mounted fix locks and our range of frames allows you to configure complex and individualized seating solutions for a wide range of children.

The dynamic back has been documented to be associated with a better postural control and less involuntary arm movements during the extensor thrusts compared to a rigid seating system (Cimolin 2009).

Read more at R82.co.uk Seating | x:panda





## **Accessories**



**Back rest extension** 8910655-x60



**Head supports** 99xxxxx-xx



Swing-away side supports 839xx-xx



Fixed side supports 8910401-x, 8910401-xP



**Upper arm supports** 8910407-x0



Overarm supports



**Vest and belts** 89113-x



Non-dynamic gas spring 891100x



Fixed rod 891150x



Quick locks 8910702-x0



**Back width extensions** 8910409-45 For x:panda size 4



Cover plate 8910406-4 For x:panda size 4

For more accessories - visit our website





Foam back plate 8910406-x



Lumbar support 8910407



**Hip supports** 8910404-4x



Armrest/Daily adjustment kit 891061x-x/8910615-11



**Extra padding** 8910404-PADx



Armrest for tray 8910730



**Grip handle for tray** 8910750



**Tray/Toy bar** 891072x-x/8910740



**Pommel** 8910408-x




**Seat wedges** 8910400-x 8910400-xx



Girth for seat cushion without abduction 891925x




**Canopy** 87401x0 Requires 9589100



## Panda Futura®

The Panda Futura Standard seat is primarily for children who need support around the arms and shoulders. This active seat is for active children who need more freedom of movement. Panda Futura is available in 5 sizes.


## **Highlights about Panda Futura**



A modular seating system with comfortable, six part cushion set



The Panda Futura seat can be mounted on many different frames. Also, the seat on a indoor frame is easy to dismount and place onto an outdoor frame.



**Back recline** 

The back can be adjusted independently of the seat. The pivotal point of recline is anatomically correct, i.e. in line with the hip joint.



25

Accessories are easily interchangeable for user comfort

Read more at R82.co.uk Seating | Panda Futura





## **Standard**

Back angling (recline), tilt in space, back and seat cushions.

## Options, **Panda Futura Active**

### Four sizes

- Size 2
- Size 2 ½
- Size 3
- Size 4

### Cushions

- Narrow
- Wide

## Options, **Panda Futura Standard**

### Five sizes

- Size 1
- Size 2
- Size 2 ½
- Size 3
- Size 4

### Cushions

- Narrow - Wide

### **Cushion colour**





### **Shell colour**







Panda Futura standard



## **Measurements**

Panda A	ctive	2	21/2	3	4
	Effective seat width	220 mm	250 mm	290 mm	320 mm
	Effective seat depth	240 mm	280 mm	310 mm	360 mm
<u>II</u>	Backrest height	350 mm	380 mm	400 mm	450 mm
	Max. load	30 kg	50 kg	50 kg	60 kg
3	Backrest angle	-4° - 27°	-10° - 20°	-9° - 27°	-6° - 17°
	Total mass	5,2 kg	5,6 kg	6,2 kg	7,5 kg
				Further meas	surements at our website

Panda St	Panda Standard		2	<b>2</b> ½	3	4
	Effective seat width	220 mm	220 mm	250 mm	290 mm	320 mm
<u> </u>	Effective seat depth	220 mm	240 mm	280 mm	310 mm	360 mm
<u>\I</u>	Backrest height	290 mm	390 mm	410 mm	440 mm	535 mm
	Max. load	30 kg	30 kg	50 kg	50 kg	60 kg
	Max. load at transportation		30 kg	50 kg	50 kg	57 kg
9	Backrest angle	-3° - 28°	-15° - 32°	-10° - 35°	-5° - 36°	-5° - 30°
	Total mass	5,7 kg	5,7 kg	6,2 kg	7,3 kg	8,7 kg
Further measurements at our website						

## Nice to know

The design of the seat and accessories ensure that the seat is easy to use and clean. The two versions; the standard and active ensure the possibility to customize the product to the activity level of the individual user.

The Panda Futura seat is popular in settings, where reuse, financially viable and less complex seating systems are warranted.

Read more at R82.co.uk Seating | Panda Futura





## **Accessories**



**Head supports** 99xxxxx-xx



**Back extension** 89655-xxx



**Back extension, curved** 89658-xxx



Wedge cushions 89760-x



Swing-away side supports



**Fixed side supports** 



Toy bar 89200





Seat extension, curved 89742-xxx



**Activity centre** 89250



**Push brace** 8997x



Push handles 89545 Panda Futura Standard 89548 Panda Futura Active

For more accessories - visit our website









Arm rest 8961x





Read more at R82.co.uk Seating | Panda Futura



## Panda Futura 5<sup>®</sup>

The Panda Futura size 5 is a comfortable supportive seat for larger children and teenagers. The seat is available in two versions; Standard and Active which provides the required comfort, support and freedom of movement.

## **Highlights about Panda Futura 5**



A range of accessories can be added



**Easy change - active to standard**Mount separate sideplates and cushions and the seat is easily transformed into the standard model. Sideplates and side cushions are accessories.



**Comfortable cushions** 



Angle adjustable back

Activate the handle on the back to carry out a smooth and stepless angle adjustment of the back.

31

Read more at R82.co.uk Seating | Panda Futura 5





## **Standard**

Available as Standard for more comfort and Active for more freedom of movement in the trunk. Angle adjustable back from -5° to 30°.

## **Options**

#### Size

- One size (Size 5)

### Cover





Gemma

Mesh







Panda Futura 5, Active

32 Seating | Panda Futura 5 R82 | a part of Etac



## **Measurements**

		Active	Standard
	Effective seat width	350 mm	350 mm
<u> </u>	Effective seat depth	430 mm	430 mm
1	Backrest height	590 mm	590 mm
	Max. load	85 kg	85 kg
	Max. load at transportation	85 kg	85 kg
9	Backrest angle	-5° - 30°	-5° - 30°
<u>_</u>	Length	560 mm	560 mm
	Overall width	560 mm	560 mm
	Total mass	15 kg	15 kg
		Further m	easurements at our website

## Nice to know

Panda Futura 5 is our seat for the largest users. It comes in two versions (standard and active) with a max. load of 85 kg and a user height of 1700 mm, which makes the Panda Futura 5 suitable for most youngsters with cerebral palsy who need adaptive seating. According to the clinical growth charts for children with cerebral palsy less than 10% of the children and youngsters weigh more than 85 kg (Brooks 2011).

Read more at R82.co.uk Seating | Panda Futura 5





## **Accessories**



**Head support** 99xxxxx-xx



**Back extension** 89658-xxxx



**Tray** 8910725-5



Back sides 89811-5xx



85430-x



**Seat extension** 



Fixed side support 8910401-4



Swing-away side support 83xxx-xxx



Arm rest 89619-x




Pommel 89930-5



For more accessories - visit our website



The pivotal point of recline is anatomically correct, being in line with the hip joint. The back angles from -5° to 30°




Seating | Panda Futura 5 Read more at R82.co.uk



# **Scallop**<sup>®</sup>

The Scallop is designed to aid floor sitting and to add extra support whilst sitting on a chair. Used by infants, children, teenagers and adults with mild to moderate disabilities as well as those without disabilities. It comes in four sizes.

## **Highlights about Scallop**



#### Removable cushion

The Scallop cushions come in two reversible colours and can be easily removed from the shells and machine-washed at  $60\,^\circ$ .



Simple and flexible



#### Portable and lightweight

Weighing less than 2kg, the Scallop is easy to travel with. You can carry the product using the two handles provided.



#### **Versatile**

The Scallop is not only used for long sitting. It can also be attached to a wide variety of chairs to add extra support while sitting around a table to enable participation in everyday activities.

37

Read more at R82.co.uk Seating | Scallop




Carry straps, seat pit, 3-point belt for size 1 and 2, 2-point belt for size 3 or 4.

# **Options**

## Four sizes

- Size 1
- Size 2
- Size 3
- Size 4

# **Cushion colours**





Red/Blue reversible cushion

Grey/Green revesible cushion

# Nice to know

Scallop provides the child with a sense of safety and security through tactile input of the seat on the body. It can be used on the floor or strapped to a chair to enhance social interaction and communication with peers.



Standard Scallop Size 1 or 2. Available in red/blue or grey/green reversible cushion



Standard Scallop Size 3 or 4. Available in red/blue or grey/green reversible cushion

38 Seating | Scallop R82 | a part of Etac



# **Measurements**

		1	2	3	4
<del>å</del>	Hip width	150-210 mm	210-240 mm	240-270 mm	270-380 mm
16	Back height	330-380 mm	380-420 mm	420-470 mm	470 + mm
	Max. load	15 kg	22 kg	29 kg	70 kg
	Total mass	0,60 kg	0,78 kg	1kg	1,44 kg
				Further meas	urements at our websit



Read more at R82.co.uk Seating | Scallop



# Wombat Living®

The Wombat Living is a cool activity chair for both toddler, pre-school, school children, teens and youngsters. The fresh look with modern colours makes it suitable for use in any indoor environment.

Available in 3 sizes.

## **Highlights about Wombat Living**



Tilt in space

The seat unit tilts from -20  $^{\circ}$  to 15  $^{\circ}$  to accommodate both relaxed and active sitting postures.



Angle adjustable push brace



**Small details - big difference** 

The cushions are made from thick comfortable foam with breathable covers. Small details like shaped cushion that protect the user's knees and extra durable edges, make a big difference in everyday use.



#### **Built-in growth**

Wombat Living offers stepless adjustment for growth. Move the backrest forward or backward and adjust the hip supports to fit any user size. The hip supports can also be angled to support the user's hips and thighs.

41

Read more at R82.co.uk Seating | Wombat Living





Tilt in space, back recline and prone angle, adjustable seat depth, adjustable back height, height and angle adjustable foot support.

# **Options**

#### Sizes

- Size 1
- Size 2
- Size 3

#### Backs

- Long back
- Short back

#### Castors

- 75 mm with brake
- 100 mm with brake
- Central brake for 100 mm castors

#### Height adjustment

- Pedal on the frame
- Release handle on the push brace
- Power actuator

### Seat angle adjustment

- Release handle beneath the seat
- Release handle on the push brace

#### Frame colour





Anthrac grey

Vhite

#### **Cushion colour**







Grey



Petroleum

contoure



Wombat Living with power actuator



Wombat Living with release handle on push brace



### **Measurements**

		1	2	3
	Effective seat width	300 mm	370 mm	420 mm
\ <u>-</u>	Effective seat depth	180-310 mm	250-380 mm	320-500 mm
<u>\I</u>	Backrest height	280-340 mm	310-380 mm	390-475 mm
<u>II</u>	Backrest height (High Back)	360-430 mm	395-460 mm	460-550 mm
2-1	Seat surface height at front edge <sup>*</sup>	360-650 mm	360-650 mm	390-670 mm
	Max. load	45 kg	45 kg	70 kg
3	Backrest angle	0°-28°	0°-28°	0°-33°
Ž.	Seat plane angle	-15°-20°	-15°-20°	-10°-22°
<u></u>	Length	600 mm	600 mm	720 mm
<u>=</u>	Overall width	550 mm	550 mm	610 mm
	Total mass**	15,5 kg	15,5 kg	22 kg

### Nice to know

Use the Wombat Living for daily activities such as reading and learning at school, having fun with friends and when helping with daily chores at home.

The height adjustable frame and the anterior tilt function ease standing transfer into the chair and ensures that unnecessary and repetitive lifting is avoided.

Although the Wombat Living at first glance looks simple and for children with mild to moderate impairments, it is possible to use the large range of accessories to adapt the chair to users with more comprehensive needs for support for the feet, knees, thighs, pelvis, trunk and head.

Read more at R82.co.uk Seating | Wombat Living





### **Accessories**



Head supports 9913092-83



Lumbar wedge 955030-x



Side supports, fixed 95121



Hip supports 95110-1, 95111-1 **Covers for hip supports** 95110-60, 95111-60



95112-1, 95113-1, 95114, 95115



Tray, smoke- coloured



Hip-belt 89661-4P



**Pommel** 



Knee abduction pommel 95230-10



Power kit 95550, 95550-3



Anti-tip device 95953



Incontinence seat cover 95562x-282

For more accessories - visit our website





Lower back extension 95140-60



Pelvis support without sides 95140-xx



**Pelvis support with sides** 95140-xx



**Back rest and seat cushion** contoured

Back rest cushion: 95503x-238 Seat cushion: 95502x-238



Adjustable push brace 95147



Helper's handles 95145



Adjustable push brace, seat mounted 95148-1, 95148-3



**Back rest and seat cushion** Back rest cushion: 95503x-2xx Seat cushion: 95502x-2xx

45





Read more at R82.co.uk Seating | Wombat Living



# **Cricket**®

The Cricket is a lightweight foldable buggy designed to suit you and your child's needs and lifestyle. This user-friendly buggy comes in two sizes and it adapts to the child's growth in several areas maximizing the number of years it can be used.

## **Highlights about Cricket**



#### **Back recline**

For extra comfort, Cricket has a two-position back recline option for the user to change from a sitting position to a relaxed position. The backrest angle can be adjusted from  $110^{\circ}$  to  $125^{\circ}$ .



**Growth adaptable** 



**Easy Fold - Light weight** 

This lightweight buggy is extremely compact when folded and requires minimum storage space.



Height & angle adjustable foot support

The Cricket comes with a unique foot support that can be height and angle adjusted to support the individual needs.

Read more at R82.co.uk Seating | Cricket 47





Back recline, adjustable seat depth, height and angle adjustable foot support, angle adjustable push brace, front wheel with swivel, rear wheel with parking brake, compact and foldable.

# **Options**

Seat

- Size 1

- Size 2

#### **Cushion colour**















Orange





Cricket standard

48 Seating | Cricket R82 | a part of Etac



### **Measurements**

		1	2
	Effective seat width	340 mm	395 mm
<u>-</u>	Effective seat depth	210-310 mm	310-430 mm
<u>\</u>	Backrest height	500 mm	570 mm
3	Seat surface height at front edge	500-530 mm	540-580 mm
*	Max. load	35 kg	50 kg
	Max. load at transportation	27 kg	43 kg
2	Backrest angle	110-125°	110-125°
>	Seat plane angle	20°	20°
>	Length	750 mm	750 mm
	Overall width	525 mm	575 mm
0	Folded length	760, mm	760 mm
M	Folded width	355 mm	355 mm
1000	Folded height	485 mm	485 mm
<b>A</b>	Total mass	12,4 kg	13,4 kg
		Further m	easurements at our website

## Nice to know

Cricket's lightweight frame is one of the most compact on the market. The frame is designed to support and ensure stability for the buggy.

The height of the push brace on the Cricket can easily be changed using the angle adjustment and the 360° turnable front wheels ensure a smooth ride and makes it easy to steer. The puncture-proof wheels provide exceptional performance. Explore the experience of the smoothest ride - from city to tougher trails, from sun to rain.

Furthermore, the Cricket is approved for transportation and complies with ISO 7176-19, for safe transportation in vehicles.

Read more at R82.co.uk Seating | Cricket





# **Accessories**



**Back rest height extension** 87043160-1, 87043160-2



**Head rest cushion** 87043150



**Sculpted head rest cushion** 87043151



**Folding top** 8704317-x, 8704327-x



**Arm rests** 87043120-1, 87043120-2



**Tray** 87043180-1, 87043180-2



**Grip bar in arm rests** 87043181-1. 87043181-2



**Grip bar in frame** 87043182-1, 87043182-2



Heel straps 87043185



Rucksack 87043183-1



**Shopping net** 87043184-1, 87043184-2



**Calf support strap** 87043190-1, 87043190-2

For more accessories - visit our website





Wedges for shaping and width reduction 87043130-x, 87043140-x



Wedges for shaping seat and back 87043140-x



Transport fittings 87043188



Directional locking mechanism, front wheels 87043187



Additional locking mechanism for foot plate 87043186



**Rain cover** 87467-210, 87467-310



Neoprene vest 89113-x



**Y-hip strap upholstered** 85427-x



**Belt harness neoprene** 89114-x



Read more at R82.co.uk Seating | Cricket



# **Stingray**<sup>®</sup>

Stingray is a superbly stylish buggy with tool free adjustment of postural support and unique swivel action for forward/rearward facing options.

The Stingray is available in 2 sizes.

### **Highlights about Stingray**



#### **Foldable**

The Stingray is folded quickly and fits into almost any car. Fold the seat and frame units seperately or altogether and pack it into your car - soon you will be on your way.



180° turnable seat unit



#### Tilt in space and back recline

Perfect seating position due to adjustable seat tilt and back angle adjustment. The seat unit tilts 45° and the back reclines down, allowing your child to rest while out and about.



#### Safe transportation in vehicles

The Stingray complies with ISO 7176-19, for safe transportation in motor vehicles. The Stingray is approved for transport in vehicles of a user weight of up to 40 kgs.

Read more at R82.co.uk Seating | Stingray 53





Tilt in space, back recline, adjustable seat depth, adjustable back height, height and angle adjustable foot support, turnable seat, frame one size.

# **Options**

#### Seat

- Size 1
- Size 2

#### Front wheels

- 12" solid or air, no swivel
- 7" solid or air, turnable front wheels, with directional stabilizer

#### Foot brake

- Individual
- Central

#### Hand brake

- Hand brakes
- No hand brakes

#### Foot rest

Long barShort bar

#### **Cushion colours**

Red







Green











Yellow

Orange







Standard Stingray with 7" turnable front wheels, with directional stabilizer

54 Seating | Stingray R82 | a part of Etac



### **Measurements**

			1	2
	Effective seat width		300 mm	350 mm
<u></u>	Effective seat depth		180-300 mm	250-400 mm
<u>\</u>	Backrest height		410-550 mm	500-650 mm
31	Seat surface height at front edge**		520 mm	520 mm
	Max. load		40 kg	40 kg
	Max. load at transportation		40 kg	40 kg
2	Backrest angle		170°	170°
7	Seat plane angle		45°	45°
<u>&gt;</u>	Length*		740 mm	740 mm
	Overall width		590 mm	590 mm
0=0	Folded length		900 mm	900 mm
M	Folded width		590 mm	590 mm
1000	Folded height		490 mm	490 mm
<u>}</u>	Total mass		19 kg	20 kg
	•	Note: The configuration of the product ** Meas	ured to adaptorplate Further m	easurements at our website

# Nice to know

The big range of accessories, flexible adjustment of the foot support and hip and side supports on the Stingray allows to make a comfortable and supportive seating solution for the individual user.

The Stingray has supreme driving properties, fantastic suspension and with the safe and stable brakes, security always has pride of place.

Read more at R82.co.uk Seating | Stingray





# **Accessories**



**Folding top** 8740110, 8740120



**Head rest cushion** 8740810, 8740820



**Sculpted head support** 872083x



Anatomic head rest 8720862, 8720863



Two-part head support




Grip handle



**Hip supports** 87402x0



Pommel 872050x



**Tray** 87200x0 8776350



**Transport fittings** 8720902




Side supports 872030x Covers for knobs on side supports 8720360



Foam back plate 877632x

For more accessories - visit our website









**Oxygen bottle holder** 8776295



Bunting bag 874013x



Vent tray 8776290



Shopping basket 8720600



Head rest fitting 8720850




Rain cover 872014x



Backpack 87207x1







# Cruiser

The Convaid Cruiser is a 30° fixed-tilt pediatric wheelchair that is convenient, lightweight and compact folding. The Cruiser is custom-built to meet each individual requirements. It comes equipped with Convaid's proprietary self-tensioning seating system, Self-Tensio™, which provides increased seating tolerance and can help promote neutral pelvic positioning.

## **Highlights about Cruiser**



**Customisable growth capabilitites**Adjustable seat depth extends chair life for a growing child.



Increased tilt angle improves postural control through upper body positioning.



#### **Foldable**

The Cruiser is folded quickly and fits into almost any car.



**Self-Tensio**<sup>™</sup>

30° fixed tilt

Built in self-tensioning seating system increases sitting tolerance.

59

Read more at R82.co.uk Seating | Cruiser





# Standard, Cruiser Standard

Solid and foldable frame, fixed tilt, anatomic back, flip-up armrest, 2-point positioning belt, off road tires, castor lock, toggle wheel locks, attendant brake and swing-away foot plates.

# Standard, Cruiser Scout

Solid and foldable frame, 30° fixed tilt, seat depth growth, one-piece angle adjustable push brace (CX12-CX16) or two-piece push handle (CX18) 2-point positioning belt, quick release rear wheels, castor lock, toggle wheel locks and removable adjustable swingaway foot plates.

# **Options**

-		
- 5	7	е

- CX10

- CX12

- CX14

- CX16 - CX18

#### Model

- Standard with foot brake
- Standard with attendant hand brake

#### Tyres

- Solid
- Solid w/transp. fitting
- Air w/transp. fitting

#### **Cushion colours**













Purple

Frame colour



Silver

Royal blue Cherry red

Panther black Pink

Green







Cruiser standard 18"



### **Measurements**

		CX 10	CX 12	CX 14	CX 16	CX 18
	Effective seat width	255 mm	305 mm	355 mm	405 mm	455 mm
_	Effective seat depth	152-299 mm	200-330 mm	254-380 mm	279-406 mm	356-533 mm
<u>II</u>	Backrest height	530 mm	560 mm	640 mm	710 mm	640 mm
<u>}</u>	Total mass	12 kg	12 kg	13 kg	14kg	14 kg
X	Max. load regular/transit	34/30 kg	34/30 kg	45,5/45,5 kg	77/77 kg	114/77 kg
	Overall width	560 mm	600 mm	650 mm	700 mm	724 mm
<u>&gt;</u>	Length	1220 mm	1120 mm	1245 mm	1320 mm	1400 mm
13	Overall height	940 mm	970 mm	940 mm	1020 mm	970 mm
0-0	Folded length	1190 mm	1220 mm	1190 mm	1350 mm	1270 mm
	Folded width	380 mm	380 mm	410 mm	410 mm	380 mm
1000	Folded height	430 mm	430 mm	430 mm	430 mm	430 mm
					Further meas	surements at our websi

# Nice to know

- The Cruiser is lightweight and has an easy three dimensional fold, which mean it requires minimal storage space and is easy to take along.
- The individual angle adjustable footplates allows for optimum leg and foot positioning.

Read more at R82.co.uk Seating | Cruiser



### **Accessories**



#### **Extended headrest cover**

With or without windows (Canopy) C904736/C904720 CX10 CX12 C904737/C904721 CX14 C904738/C904722 CX16 C904739/C904723 CX18 C904740/C904724



Rain cover

CX10-14 C904939 CX16-18 C904910



#### **Headrest extension**

Curved cordura (black) CX10 C903332 CX12 C902863 CX14 C900484 CX16 C900138 CX18 C902619



#### Headwings

Padded: CX10 C903331 CX12 C902862 CX14 C900681 CX16 C900142 CX18 C902618 Occi: C904779



3-pt positioning belt

With depth adjustement croth strap



H-harness with padded covers

3-pt position belt required



Full support torso vest

CX10 C903330 CX12 C902861 CX14 C900680 CX16 C900339 CX18 C902616

11



Soft lateral support

Single/Double flap C903328/C903329 CX12 C902859/C902860 CX14 C900678/C900679 C900140/C900204 CX16 C902614/C902615 CX18



Reducer seat insert

Seat width 5 cm smaller, seat depth 2,5

cm smaller

CX10 C904211 C904111 CX12 CX14-18 C904110



**Calf panel** 

C904580



Medial thigh support-abductor

with position/align cushion CX10 C903326 CX12 C902857 C900480 CX14 CX16 C900137 CX18 C902607



Lateral thigh support-abductor

with position/align cushion CX10 C903327 CX12 C902858 C900482 CX14 C900188

CX18 C902608





**Angle adjustable footplates**CX10 C903445
CX12-14 C903268

CX16-18 C903269



Footplate securement strap C903584



Foot positioners

Small (18-25 cm) C903415 Medium (20-28 cm) C903416 Large (23-30 cm) C903417



Upper extremity support surface

Clear tray CX10 C903078 CX12 C903385 CX14 C903386 CX16 C903387 C903388 CX18

16




Anchor transit kit

CX10-16 C997125 CX18 C997141 Scout model: CX12-16 C997219 CX18 C997220



Anti-tip

CX10-12 C904520 CX14-16 C904540



Under seat storage basket

CX10-12 C904520 CX14-16 C904540 CX18 C904048



**Utility** bag C904531



Saddle bag C904094



Wheelchair travel bag

C902022



**Curb tipper** 


CX10-16 C904128 CX18 C903916 C903916 Scout:



Incontinence cushion liner

63

C902324



# **EZ Rider**

The Convaid EZ Rider is an upright, 10° fixed tilt pediatric wheelchair that is best-in-class, convenient, lightweight and compact folding. The EZ Rider is custom-built to meet each individual requirements.

## **Highlights about EZ Rider**



**Customisable growth capabilitites**Adjustable seat depth extends chair life for a growing child.



**10° fixed tilt**Provides upright posture to enhance functional positioning.



**Upholstery & colours**Removable, washable Cordura available in various colours.



**Foldable**The EZ Rider is folded quickly and fits into almost any car.

Read more at R82.co.uk Seating | EZ Rider 65





 $Removable\ adjustable\ swing-away\ footplates,\ one-piece,\ angle\ adj.\ push\ handle\ (CX12/14/16),\ 2-Piece\ Push\ Handles\ (CX18),$ 30-degree fixed tilt, 2-point positioning belt, seat depth growth of 13 cm, individual footbrake.

# **Options**

- EZ14 - EZ16

- EZ18

### - Solid

- Air

- Solid w/transp. fitting

- Air w/transp. fitting

#### **Cushion colours**















Frame colour



Silver







EZ Rider Standard, size 18"



### **Measurements**

		EZ 12	EZ 14	EZ 16	EZ 18
	Effective seat width	305 mm	355 mm	405 mm	455 mm
<u>_</u>	Effective seat depth	230-360 mm	280-410 mm	305-430 mm	370-495 mm
1	Backrest height	530 mm	610 mm	660 mm	680 mm
2	Total mass	11 kg	14 kg	14kg	13 kg
	Max. load regular/transit	34/30 kg	45,5/45,5 kg	77/77 kg	114/77 kg
	Overall width	546 mm	548 mm	660 mm	660 mm
<u>&gt;</u>	Length	915 mm	950 mm	1105 mm	1220 mm
1	Overall height	965 mm	980 mm	1015 mm	950 mm
0-0	Folded length	1130 mm	1170 mm	1220 mm	1160 mm
	Folded width	360 mm	340 mm	360 mm	370 mm
1000	Folded height	445 mm	460 mm	430 mm	390 mm
				Further m	easurements at our websit

# Nice to know

 EZ Rider is Lightweight, portable and folds compactly for the active user. It is suitable for children in need of a functional upright seating position, where the focus is on allowing the child to active and aware of the surroundings.

Read more at R82.co.uk Seating | EZ Rider



### **Accessories**



**Extended headrest cover** 

With or without windows (Canopy) C904737/ C904721 EZ12 EZ14 C904741/ C904725 EZ16 C904742/ C904726 EZ18 C904743/ C904727



**Headrest cover** 

C900853 EZ12 EZ14 C901210 EZ16 C901205 EZ18 C901211



Rain cover and mosquito net

Mosquito net: 12/14 C904655 Rain cover: 12/14 C904939 16/18 C904940 16/18 C904656



**Headrest extension** 

Curved cordura (black) EZ12 C902863 EZ14 C900484 EZ16 C900138 EZ18 C902619



Headwings

Padded: C900849 EZ14 C900682 FZ16 C900683 C900684 EZ18

Occi: C904779



H-harness with padded covers

3-pt position belt required



3-pt positioning belt

With depth adjustement croth strap



Full support torso vest

EZ12 C900848 EZ14 C900338 FZ16 C901020 FZ18 C900387



Soft lateral support

Single/Double flap

C900846/ C900847 EZ12 C900327/ C900650 EZ14 C901018/ C900651 EZ16 C900386/ C900652



**Calf panel** 

C904580



Reducer seat insert

Seat width 5 cm smaller, seat depth 2,5 cm smaller

EZ12 C904211 EZ14-18 C904110



Medial thigh support-abductor

with position/align cushion

EZ12 C900843 C900329 EZ14 EZ16 C901024 C900388 EZ18





Lateral thigh support-abductor

with position/align cushion EZ12 C900844 EZ14 C900330 C901026

C900389

EZ16

EZ18



**Angle adjustable footplates**EZ12/14 C903268
EZ16/18 C903269 EZ12/14 EZ16/18



Footplate securement strap

C903584

**1**5



**Foot positioners** 

Small (18-25 cm) C903415 Medium (20-28 cm) C903416 Large (23-30 cm) C903417



Height adjustable

**flip-up armrests**EZ12-14 C904523
EZ16-18 C904524



Upper extremity support surface

Clear tray EZ12 C904525 C904526 EZ14 EZ16 EZ18 C904527 C904528



Upper extremity support surface

Foam tray EZ12 C903669 EZ14-16 C903670 EZ18 C903671



69

**Anchor transit kit** C997267



Anti-tip EZ12-14 C900377 EZ16-18 C900378



Utility bag C904531



Wheelchair travel bag C902022

Read more at R82.co.uk Seating | EZ Rider



# Rodeo

The Convaid Rodeo is a popular pediatric tilt-in-space wheelchair that is convenient, lightweight and folds compactly. The Rodeo provides up to 45° of adjustable tilt that offers pressure relief and upper body support. This tilt-in-space Rodeo is fully customizable to meet user's specific therapeutic needs and provides growth option through its adjustable seat depth feature.

### **Highlights about Rodeo**



Adjustable tilt and back recline

Provides pressure relief and support.



**Customisable growth capabilitites**Adjustable seat depth extends chair life for a growing child.



Lighweight and durable

Lightweight frame combined with high Max. load - up to 77 kg.



**Individual positioning** 

Accommodates functional or rest positions and increased sitting tolerance for maximum comfort.

71

Read more at R82.co.uk Seating | Rodeo



5 to 45 degrees of adjustable tilt, recline, central mounted foot brake, adjustable one-piece push handle, closure strap, seat depth growth of 13 cm, frontwheels with directional stabilizer, 2-point positioning belt, elevating legrest /one-Piece, angle adjustable footplate.

# **Options**

#### Size:

- RD12
- RD14
- RD16

#### Brake

- Foot brake
- Attendant hand

#### Tyres

- Solid - Air
- Solid
- w/transp. fitting
- Air w/transp. fitting

#### **Cushion colours**

















Frame colour



Silver







Folded Rodeo



# **Measurements**

		RD 12	RD 14	RD 16
	Effective seat width	305 mm	355 mm	405 mm
<u>_</u>	Effective seat depth	230-360 mm	280-410 mm	330-460 mm
<u>II</u>	Backrest height	560 mm	610 mm	635 mm
	Total mass	15 kg	16 kg	17 kg
	Max. load regular/transit	34/30 kg	46/46 kg	77/64 kg
	Overall width	580 mm	630 mm	680 mm
7	Length	890 mm	890 mm	910 mm
13	Overall height	1090 mm	1120 mm	1140 mm
0	Folded length	860 mm	860 mm	1015 mm
	Folded width	585 mm	635 mm	685 mm
1000	Folded height	430 mm	460 mm	460 mm
			Further m	easurements at our website

# Nice to know

- Rodeo is suitable for walks in uneven terrain and nature
- It is easy to adjust most of the settings and it is possible to use it for different children.

Read more at R82.co.uk Seating | Rodeo



### **Accessories**



**Extended Headrest cover** 

With or without windows (Canopy) 12 C904745/ C904729 14 C904746/ C904730 16 C904747/ C904731



Rain cover and Mosquito net

Mosquito net: 12/14 C904655 Rain cover: 12/14 C904939 C904940 16 C904656



**Headrest Extension** 

Required for transit model 12 C902863 14 C900484 C900138



Headwings

Padded: C902673 14 C902653 16 C903147 Occi: C904779



H-harness with padded covers

With padded covers C903988



3-pt positioning belt

With adjustable crotch strap

C904017



Full torso vest

C903261 14 C903260 16 C903259



Reducer seat insert

Seat width 5 cm smaller, seat depth 2,5

C904211 12

14 C904110



Swing away laterals - contoured

X-small C902980 Small C902981 Medium C902982 Large C902983





**Cushions** 



Calf panel C904361



Thigh support

Lateral: Medial: 12 C903668 12 C902668 14 C903648 14 C902647 C903141 16 C903140 16









Height adjustable flip-up armrest C904346





Rear anti-tip tubes (pair)



Upper extremity support surface (clear tray)

Flip-up arm rest required 12 C903405 14 C903401 16 C903402





Anchor transit kit



Medical necessity storage 12 C904476 14/16 C904477





Oxygen tank bag



Utility bag









Environment Indoor & outdoor\*

# **High-low Frame**<sup>®</sup>

The High-low Frame is a stable base with a user friendly height adjustment from floor to table height.

The High-low Frame is available either as an indoor frame with small lockable castors or as an outdoor frame with large wheels and foot brake.

# **Highlights about High-low Frame**



### **Height adjustment**

The height is adjustable from floor to table height. The frame are available with either gas or power height adjustment.



#### Interface

Various seats can be mounted on the High-low Frame by means of an interface. Here a Hich-low Frame for outdoor mounted with a Panda Futura.



### **Angling of seat base**

The seat base has an angle range of  $40^\circ$  which helps obtain varied seating positions. The base can be set to angle from  $0^\circ$  to  $30^\circ$  or from  $-10^\circ$  to  $20^\circ$ .



#### Safe transportation in vehicles

The High-low Frame, outdoor, complies with the international transit wheelchair standard ISO 7176-19. The High-low Frame, outdoor, can be used with seat units that comply with ISO 16840-4.

77

Read more at R82.co.uk Seating | High-low Frame





Height adjustment, angling of seat plane, one size outdoor frame with 12  $\frac{1}{2}$ " rear wheels 7" swivel front wheels and foot brake, one size indoor with 100 mm lockable castors.

### **Options**

#### Indoor frame

- Gas
- Power

#### Outdoor frame

- Gas
- Power
- Handbrakes
- Handbrakes and power

### Frame colour



Black


### **Accessories**



Push brace for indoor High-low



Angle adjustable push brace for indoor frame 8910565, 8910555



**Extra gas spring** 89955+8911003



Anti-tip for High-low 89951-0x



Power kit for High-low 8995x



Rod for foot support Double rod: 891084x-x Single rod: 891082x-x



**Parted foot supports** 89146-x



**Foot plate** 95145-1, 85145-2

For more accessories - visit our website





Individual angle adjustable foot supports 89143-x



**Foot plate model Wombat** 95145-2, 95145-3, 95145-4



Individual multi adjustable foot support Size 1 8910851-x Size 2+3 8910852-x

foot Adj




Adjustable legsupport 89162



Foot plate 89130



Foot plate, padded 89128



Transport fittings only for outdoor frames 9809xx



Wheels with directional stabilizers 89592-KIT

79

## **Measurements**

		Indoor Frame		<b>Outdoor Frame</b>		
		w/Gas	w/Power	w/Gas	w/power	
<b>.</b> Z. <b>1</b>	Seat surface height at front edge*	170-600 mm	260-700 mm	260-700 mm	260-700 mm	
<u>*</u>	Max. Load	60 kg	60 kg	60 kg	60 kg	
	Max. load at transportation	-	-	43 kg	43 kg	
Ź.	Seat plane angle	-10° - 30°	-10° - 30°	-10° - 30°	-10° - 30°	
Σ.	Length	690 mm	690 mm	960 mm	960 mm	
≞	Overall width	560 mm	560 mm	620 mm	620 mm	
<u></u>	Total mass	12,5 kg	13 kg	17,5 kg	18,0 kg	
*measured to adaptor plate Further measurements at our website						

Read more at R82.co.uk Seating | High-low Frame



# **High-low:x Frame**<sup>®</sup>

The High-low:x Frame is a stable base with a user friendly height adjustment from floor to table height.

The High-low:x Frame is available either as an indoor frame with small lockable castors or as an outdoor frame with large wheels and foot brake.

# **Highlights about High-low:x Frame**



### **Height adjustment**

The height is adjustable from floor to table height. The frames are available with either gas, power or hydraulic height adjustment.



### **Fixed mounting or interface**

On the High-low:x size 1 and 2 the seats are mounted into the interface. On size 3 the seats can be mounted either directly onto the High-low:x Frame or by using an interface.



### **Angling of seat base**

The seat base has an angle range of  $40^{\circ}$  which helps obtain varied seating positions. The base can be set to angle from  $0^{\circ}$  to  $30^{\circ}$  or from  $-10^{\circ}$  to  $20^{\circ}$ .



### Safe transportation in vehicles

The High-low:x Frame, outdoor, complies with the international transit wheelchair standard ISO 7176-19. The High-low Frame, outdoor, can be used with seat units that comply with ISO 16840-4.

81

Read more at R82.co.uk Seating | High-low:x Frame





Height adjustment, angling of seat plane, 2 sizes outdoor frame with  $12 \frac{1}{2}$ " rear wheels 7" swivel front wheels and foot brake, 3 sizes indoor with 100 mm lockable castors.

# **Options**

#### Height adjustment

- Gas
- Power
- Hydraulic

#### Seat mounting

- Adaptor plate
- Fixed mounting on seat plate (Only for size 3 indoor)

#### Frame colour



Black



Standard High-low:x Frame with gas



High-low:x Frame mounted with x:panda



## **Measurements**

High-low	High-low:x indoor		2 gas	1 power	2 power	3 power	2 hydraulic
.Z. <b>t</b>	Seat surface height at front edge*	170-600 mm	250-610 mm	170-600 mm	250-620 mm	320-560 mm	220-580 mm
<u></u>	Max. load	60 kg	60 kg	60 kg	60 kg	100 kg	60 kg
Œ.	Seat plane angle	-10°- 30°	-10°- 30°	-10°- 30°	-10°- 30°	-10°- 30°	-10°- 30°
Σ.	Length	690 mm	780 mm	690 mm	780 mm	780 mm	780 mm
盡	Overall width	560 mm	620 mm	560 mm	620 mm	620 mm	620 mm
.Z.	Total mass	12,5 kg	14,0 kg	13,0 kg	16,5 kg	17,0 kg	16,5 kg
				*To	adaptor plate Fu	urther measureme	nts at our website

High-low	:x outdoor	1 gas	2 gas	1 power lift	2 power lift	2 hydraulic
.Z. <b>t</b>	Seat surface height at front edge*	260-700 mm	330-640 mm	260-700 mm	330-640 mm	430-590 mm
<u>.</u>	Max. load	60 kg	60 kg	60 kg	60 kg	60 kg
	Max. Load at transportation	43 kg	43 kg	43 kg	43 kg	43 kg
Σ.	Seat plane angle	-10°-30°	-10°- 30°	-10°- 30°	-10°- 30°	-10°- 30°
Σ.	Length	960 mm	1110 mm	960 mm	1110 mm	1110 mm
益	Overall width	620 mm	620 mm	620 mm	620 mm	620 mm
<u>Z.</u>	Total mass	17,5 kg	20,0 kg	18,0 kg	20,5 kg	20,5 kg
				*To adaptor pla	ate Further measure	ments at our website

# Nice to know

Recommendations for the choice of frame.

Stingray seat size 1 or 2: Frame size 2 x:panda size 1-2: Frame size 1 x:panda size 3: Frame size 2 or 3 x:panda size 4: Frame size 3

User (weight of child) + Seat + accessories = Load Example: Child (25 kg), x:panda size 2 (8.5 kg), Accessories (5 kg), Total: 38.5 kg.

Read more at R82.co.uk Seating | High-low:x Frame





# **Accessories**



Push brace for indoor High-low:x 8910545



Angle adjustable push brace for indoor frame 8910565, 8910555



Foot pedal kit 8910610



Extra gas spring 89955+8911003



Power kit for High-low:x



**Anti-tip for High-low:x** 89951-0x



Transport fittings only for outdoor frames
9809xx



Double rod for foot support



Single rod for foot support 891082x-x



Parted foot supports 89146-x




**Foot plate** 95145-1, 85145-2



Individual angle adjustable foot supports 89143-x

For more accessories - visit our website





Individual multi adjustable foot **support** Size 1 8910851-x Size 2+3 8910852-x



**Adjustable legsupport** 89162



**Foot plate model Wombat** 95145-2, 95145-3, 95145-4



Foot plate 89130



Foot plate, padded 89128



Wheels with directional stabilizers 89592-KIT



Central brake kit 8910950







Max. load 100 kg

Environment Indoor & outdoor

# Multi Frame®

The Multi Frame is an innovative, lightweight and highly manoeuvrable wheelchair base designed to fit both the Panda Futura and x:panda, but also other manufacturers seating that is constructed in a similar fashion and meets the current standards.

# **Highlights about Multi Frame**



Tilt in space

Simple to operate tilt in space with the minimum of moving parts.



**Seats** 

Mount a Panda Futura or x:panda seat on the Multi Frame.



**Prepared for other manufactures seats and backs** 

Mount a Panda Futura or x:panda seat on the Multi Frame.



### Safe transportation in vehicles

The Multi Frame complies with ISO 7176-19, transporting users in motor vehicles. The Multi Frame can be used with seat units that complies with ISO 16840-4.

87

Read more at R82.co.uk Seating | Multi Frame





3 sizes, built-in tilt in space 0 to 30°, change of seat to floor height, brakes, anti tips - swing away type, angle adjustable push brace.

# **Options**

#### Three sizes

- Size 1
- Size 2 - Size 3

#### Seat mounting

- Fixed mounting on seat plate (standard)
- Adaptor plate

#### Rear wheels

- 12 ½" solid or air
- 22" solid or air

#### Brakes

- Wheel brake
- Drum brake

#### Front wheels

- 7" solid grey or black
- 7" air grey or black

### **Framecolour**



Silver



Multi Frame with 12.5" rear wheels



Multi Frame with 22" rear wheels

88 Seating | Multi Frame R82 | a part of Etac



# **Measurements**

8,	Seat surface height at front edge <sup>+</sup>	005 405 455		
. A OR		395, 425, 455 mm	395, 425, 455 mm	395, 425, 455 mm
	Max. Load	100 kg	100 kg	100 kg
O . U	Max. load at transportation	100 kg	100 kg	100 kg
3,	Seat plane angle	0-30°, 35°, 40°	0-30°, 35°, 40°	0-30°, 35°, 40°
	Width between wheels	340 mm	390 mm	440 mm
<u>P</u> ,	Length 12½"	640 mm	640 mm	690 mm
<u>P</u> ,	Length 22"	750 mm	750 mm	800 mm
<u> </u>	Overall width	590 mm	640 mm	690 mm
19,	Overall height**	915 mm	915 mm	915 mm
<u>G</u> ,	Folded length 12½"	730 mm	730 mm	780 mm
-	Folded width 12½"	470 mm	470 mm	470 mm
ig,	Folded height 12½"	640 mm	640 mm	690 mm
<u>G</u> ,	Folded length 22"	850 mm	850 mm	900 mm
<u> - </u>	Folded width 22"	590 mm	640 mm	690 mm
ig,	Folded height 22"	560 mm	560 mm	560 mm
<u>₽</u> 8₩	Total mass 12 ½" wheels	13 kg	13 kg	13 kg
<u>P</u> 8	Total mass 22" wheels	14 kg	14 kg	14 kg

Read more at R82.co.uk Seating | Multi Frame





## **Accessories**



**Interface** 8911610-98 + 8911610-99



Transport fittings 970905



Spoke protector 970922



Single rod for footsupport 8910821-x



Individual multi adjustable foot support 8910852-x



**Foot plate** 95145-1, 85145-2



**Foot plate model Wombat** 95145-2, 95145-3, 95145-4



Vent tray

For more accessories - visit our website

# Nice to know

The angle of the push brace can be adjusted to make a bigger distance between the push brace and seat. The adjustment is easily done with the use of the strap under the frame and the correct position is easily recognized by the click of the locking pin into the fitting.

90 Seating | Multi Frame R82 | a part of Etac











Overall width 660 mm

Max. load 100 kg

Environment Indoor & outdoor

# Multi Frame:x®

Designed to provide extraordinary stability and safety, the Multi Frame:x is the perfect choice of wheelchair base for large seat units like the Panda Futura 5 and x:panda 4. The Multi Frame:x is a strong, stable and highly manoeuvrable wheelchair base.

## **Highlights about Multi Frame:x**



### **Seats for Multi Frame:x**

Multi Frame:x is the perfect choice of wheelchair base for large seat units like the Panda Futura 5 and x:panda 4 from R82.



### **User-friendly adjustability**

The tilt-in-space and back recline functions are operated from the push brace ensuring optimum working conditions. Stepless recline of the push handle makes it possible to adapt the height to any carer.



#### Seat height adjustment

Positioning the user as close to the floor as possible ensures a much easier transfer.



#### **Tested and approved**

The Multi Frame:x complies with ISO 7176-19, transporting users in motor vehicles. The Multi Frame:x can be used with seat units that comply with ISO 16840-4.

93

Read more at R82.co.uk Seating | Multi Frame:x





Built-in recline, built-in tilt, anti-tip device, seat height adjustment, angle adjustable push brace.

# **Options**

#### Rear wheels

- 12 ½" solid or air
- 20"/22"/24" solid

#### Front wheels

- 6" solid
- 7" solid or air

#### Brakes

- Wheel brake
- Drum brake

### **Framecolour**



Anthracite grey



Multi Frame:x with 12  $\frac{1}{2}$ " rear wheels and 7" front wheels



Multi Frame:x with 22" rear wheels and 6" front wheels



### **Measurements**

		3
Pri	Seat surface height at front edge	410, 435, 460 mm
& <sub>s</sub>	Max. Load	100 kg
	Max. load at transportation	100 kg
<i>B</i> ,	Seat plane angle	-5-30°
	Width between wheels	215 x 405 mm
<u>P</u> ?	Length	840 mm
	Overall width	660 mm
P	Overall height	980 mm
<u>B</u> 8▲	Total mass	20 kg

The Multi Frame:x is approved for transporting users in motor vehicles. Secure the user with an approved 3-point belt according to ISO 10542-1 and use an approved 4-point strap-type tiedown system according to ISO 10542-2 to secure the wheelchair in a forward-facing position in the vehicle. All removable accessories must be dismounted and secured separately in the vehicle.

# Nice to know

Use the seat height adjustment to position the user as close to the floor as possible for easier transfer – or use the seat height extension to increase the distance from the frame to the floor to 460 mm, which allows for customisation for users with longer legs.

Read more at R82.co.uk Seating | Multi Frame:x





# **Accessories**



Adapter 8911610-99L



Mud guards 96052x



**Spoke protector** 9709xx



**Seat height extension** 961030-xx



Mounting for foot rest system



Single rod for foot rest



**Foot plate model Wombat** 95145-2, 95145-3, 95145-4



**Foot plate** 95145-1, 85145-2



**Calf support** 85131



**Double rod for foot rest** 897528x



**Parted foot supports** 89146-x



Individual angle adjustable foot supports 89143-x

For more accessories - visit our website









**Transport fittings** 960905



**Vent tray** 95281-10













Overall width 590 mm

Max. load 40 kg

Environment Indoor & outdoor

# Stingray Frame®

The Stingray frame with 180° turntable base makes it possible to move your child from forward to rearward facing without leaving their seat.

The carbon Stingray Frame is in a league of its own with the trendsetting design, great adjustments and ease of steering.

# **Highlights about Stingray Frame**



#### 180° turnable seat unit

Revolutionary function with a 180  $^\circ$  turnable seat unit. Let your child explore the surroundings during the trip or simply turn the unit to obtain eye contact with your child.



#### Foldable frame

The Stingray Frame is foldable for storage and transportation.



#### Tilt in space

The tilt in space function is very useful when the user needs to rest. Simply angle the complete unit into the desired position.



#### Safe transportation in vehicles

The Stingray Frame complies with the international transit wheelchair standard ISO 7176-19. The Stingray Frame can be used with seat units that comply with ISO 16840-4.

99

Read more at R82.co.uk Seating | Stingray Frame





180° degree turntable base, one size with 12" rear wheels, tilt-in-space, adaptor plate.

### **Options**

#### Front wheels

- 12" solid or air, no swivel
- 7" solid or air, turnable front wheels, with directional stabilizer

#### Foot brake

- Individual
- Central

#### Hand brake

- Hand brakes
- No hand brakes

#### Foot rest

- Long bar
- Short bar

#### Frame colour



Black

## **Accessories**



**Shopping basket** 8720600



Transport fittings 8720902



Vent tray 8776290



Oxygen bottle holder

### For more accessories - visit our website



Standard Stingray with 12" no swivel front wheels



Standard Stingray with 7" turnable front wheels, with directional stabilizer



# **Measurements**

		1
3	Seat surface height at front edge**	520 mm
	Max. Load	40 kg
	Max. load at transportation	40 kg
>	Seat plane angle	170°
<u>}</u>	Length*	740 mm
	Overall width	590 mm
13	Overall height	1000, 1030, 1060 mm
0-0	Folded length	900 mm
	Folded width	590 mm
0-0	Folded height	490 mm
	Total mass	20 kg
	*Measured with turnable front w	wheels ** Measured to adaptorplate Further measurements at our website





101

Standard Stingray with 7" turnable front wheels, mounted with x:panda seat

Read more at R82.co.uk Seating | Stingray Frame









Max. load 100 kg




Environment Indoor & outdoor

# Combi Frame:x®

The Combi Frame:x is a unique wheelchair frame for customized seating systems. The Combi Frame:x is ideal for both indoor and outdoor use and can be fitted with different seat and back systems. The Combi Frame:x is suitable for users needing an individual seating system and base.

# **Highlights about Combi Frame:x**



### Centre of gravity and tilt in space

The Combi Frame:x centre of gravity is self-centring. The seating unit can be tilted from horizontal to 40 degrees.



# **Compatible with different seat and back systems**

The Combi Frame:x can be equipped with different back systems and cushions.



### **Fixed mounting or interface**

Various seats can be mounted directly onto the Combi Frame:x. Alternatively the seats can be mounted into an interface. The frame is used for a variety of seating units like the Panda Futura and x:panda from R82.



#### Safe transportation in vehicles

The Combi Frame:x complies with the international transit wheelchair standard ISO 7176-19. The Combi Frame:x can be used with seat units that comply with ISO 16840-4.

Read more at R82.co.uk Seating | Combi Frame:x 103





4 sizes, self-centring centre of gravity, prepared for x:panda, Panda and moulded seats, anti-tip device.

# **Options**

#### Four sizes

- Size 1
- Size 2
- Size 3
- Size 4

#### Seat mounting

- Fixed mounting on seat plate (standard)
- Adaptor plate

#### Rear wheels

- 12 ½" solid or air
- 20"/22"/24" solid or air

#### Front wheels

- 6" solid or air
- 7" solid or air

#### Back recline

- No recline
- Recline

#### Brakes

- Wheel brake
- Drum brake
- Wheel brake + drum brake

#### Frame colour







White

lver Black



Combi Frame:x with 121/2" wheels



Combi Frame:x with 24" wheels

104 Seating | Combi Frame:x R82 | a part of Etac



### **Measurements**

		1	2	3	4
921	Seat surface height at front edge	430, 460, 490 mm			
8	Max. Load	100 kg	100 kg	100 kg	100 kg
8	Max. load at transportation	100 kg	100 kg	100 kg	100 kg
<i>B</i> '	Seat plane angle	0-40°	0-40°	0-40°	0-40°
	Width between wheels	440 mm	490 mm	560 mm	600 mm
<u>P</u>	Length	830 mm	830 mm	830 mm	830 mm
<u> </u>	Overall width	580 mm	630 mm	680 mm	730 mm
19°	Overall height	900 mm	940 mm	1040 mm	1140 mm
000	Folded length	800 mm	800 mm	840 mm	940 mm
	Folded width	580 mm	630 mm	680 mm	730 mm
10000	Folded height	550 mm	550 mm	550 mm	550 mm
<u></u> 8 ■	Total mass	22 kg	23 kg	24 kg	27 kg
				Further meas	surements at our website

# Nice to know

The Combi Frame:x is a unique wheelchair frame for customised seating systems. The multiple adjustment possibilities and the large accessories selection allows for a highly customised seating solution for the individual user with complex needs.

The kit for dynamic back rest with recline provides dynamic movement of the back rest. A dynamic back allows the user to push back during the extensor trust and supports the user to regain the starting point for the movement.

Read more at R82.co.uk Seating | Combi Frame:x





# **Accessories**



Mounting for head support



Mounting for head support 9802x0



Kit for dynamic back rest with recline 980946-DYN



Fix locks 970710



Mudguards 98052x-1



**Spoke protectors** 97102x



**Transport fittings** 980900, 980903



8911610-98 + 8911610-99



Mounting for single/parted foot support 98015x



Single rod for foot rest



Bar for parted foot support 897528x



Individual angle adjustable foot supports 89143-x

For more accessories - visit our website





Foot support 89144-x01



**Parted foot supports** 89146-x



**15** 

**Foot plate model Wombat** 95145-2, 95145-3, 95145-4



**Foot plate** 95145-1, 85145-2



Leg rest with gas strut 970500-x



Fixed angle leg rest



Calf supports, parted incl. cover



Calf support, one piece incl. cover 970555-x

21




**Extra cover for calf support** 970550-x13, parted 970555-x13, one piece



Individual, multi-adjustable foot support 981085x-x



Vent tray 970815



Push brace, complete (Supplied as standard for Combi Frame:x as of 1st January 2019) 97054x-1









Overall width 600 mm

Max. load 120 kg

Environment Indoor only

# **Strong Base**<sup>®</sup>

The Strong Base is a simple, stable and sturdy base for large seat units such as the Panda Futura 5 and the x:panda 4.

Seat units are mounted either directly onto the Strong Base or by means of an interface making it easy to switch between frames.

## **Highlights about Strong Base**



## **Height adjustment**

The Strong Base has either manual or spindle height adjustment. The manual base adjust from 30 to 55 cm in height - the spindle base adjust from 30 to 51 cm in height.



#### Tilt in space

The tilt in space function helps obtain varied seating positions. Simply angle the complete unit into the desired position. Tilt in space from -10 to  $+25\,^{\circ}$ .



## Fixed mounting or interface

Various seats can be mounted directly onto the Strong Base. Alternatively the seats can be mounted by means of an interface.



#### **Easy to manoeuvre**

The smooth running castors make it very easy to move the Strong Base around. All four castors are equipped with directional stabilizers and brakes.

Read more at R82.co.uk Seating | Strong Base 109





1 size in 2 versions, built-in tilt in space, 125 mm swivel castors with brakes and directional stabilizers, lock on all 4 wheels, seat plate.

## **Options**

#### Seat mounting

- Seat plate (standard)
- Adapter plate

#### Height adjustment

- Manual
- Spindle

#### Frame colour



Anthracite grey

## **Accessories**



**Push brace** 89105x5



Angle adjustable push brace for indoor frame



Adapter for seat mounting



Single rod for foot supports



**Double rod for foot suppports** 891084x-x



Individual multi adjustable foot Size 1 8910851-x Size 2+3 8910852-x



Individual angle adjustable foot supports 89143-x



Wombat foot plate 95145-3, 95145-4

For more accessories - visit our website





**Parted foot supports** 89146-x

Width extension for multiadjustable foot support 8910855

#### **Measurements**

		1	2
921	Seat surface height at front edge <sup>+</sup>	420-670 mm	420-630 mm
8	Max. Load	120 kg	120 kg
<b>E</b> ,	Seat plane angle	-10°-25°	-10°-25°
<u>P</u> ?	Length	700 mm	700 mm
<u> </u>	Overall width	600 mm	600 mm
192	Overall height	300-550 mm	300-510 mm
<u>888</u>	Total mass	14 kg	15,5 kg
		*Measured from x:panda 4 seat to	floor Further measurements at our website

## Nice to know

The Panda Futura 5 or x:panda size 4 is mounted directly onto the Strong Base or by means of an interface making it possible to switch the seat between frames.

The seat surface height at the front edge of the x:panda size 4 on 420 mm allows standing transfers into the seat.

Read more at R82.co.uk Seating | Strong Base





# **Standing, overview**

		Meerkat	Rabbit Up	Toucan	Caribou	Gazelle
Target users	Max. load	90 kg	70 kg	70 kg	100 kg	80 kg
	User height	700-1800 mm	1000-2000 mm	1000-2000 mm	850-1840 mm	700-1350 mm
Features	Easy transfer using the tilt function				•	•
	Accessories are easily interchangeable	•			•	
	Mobile (Quick release wheels)		•			
	Hip abduction (0-30°)					•
	Dynamic standing	•				
Environment	Indoor	•	•	•	•	•
	Outdoor		•			

112 R82 | a part of Etac



# X Standing

Safe support while standing increases both physical and mental well-being!



# **Meerkat**®

The Meerkat is a simple and modular built standing frame with a free choice of positioning in anterior or posterior.

The Meerkat is available in 3 sizes and suitable for children and youngsters.

## **Highlights about Meerkat**



#### **Dynamic standing**

The child's freedom of movement is encouraged with minimal support, pure product design and the innovative rocker.



#### Simple and intuitive

All support is intuitively easy to use and adjust, making life easier for parents and carers. The low entrance height ensures easy transfer.



#### **Easy adjustments**

Simple adjustments can be done without any tools required. The child can have fun and play while standing.



#### **Growth adjustable accessories**

Growth opportunity and increased reuse possibility.

Read more at R82.co.uk Standing | Meerkat 115





Base plate with wheels, column, tool.

## **Options**

Size

- Size 1
- Size 2
- Size 3


## **Base plate colours and castors**



Strawberry ice with cow, size 1




Pastel blue with airplane,



Anthracite with grey, size 2+3







116 Standing | Meerkat R82 | a part of Etac



## **Measurements**

		1	2	3
<b>Y</b> I	User height: floor to axil	400-750 mm	720-1000 mm	840-1300 mm
	User height	700-1050 mm	950-1350 mm	1170-1800 mm
	Max. load	20 kg	65 kg	90 kg
•	Length	595 mm	800 mm	800 mm
<u></u>	Overall width	525 mm	780 mm	780 mm
[ <u>.</u>	Overall height	375-750 mm	500-1000 mm	650-1300 mm
<u>.å.</u>	Total mass	9 kg	15 kg	15 kg
				Further measurements at our websi

## Nice to know

The Meerkat can be used to create a learning environment, where the chest, pelvic and knee supports are used to provide the child with the exact support needed to challenge static and active trunk control. When the Meerkat is fitted with a rocker, it is possible to challenge the child's reactive trunk control.

When the child develops control over a segment of the trunk, the support can gradually be moved from top down.

Thereby the Meerkat allows the child to develop trunk control, as seen in typically developing children, which might lead to improvement in gross motor function and mobility.

Read more at R82.co.uk Standing | Meerkat





## **Accessories**



Chest support 86948x



Pelvic support 86924x



Knee supports 86984x



Sandals 86949x-x



**Rocker** 118598-1, 117047-1



**Tray** 869485



Displacement kit



Daily adjustment kit 124026



For more accessories - visit our website

118 Standing | Meerkat R82 | a part of Etac



Easy & fun

Mounting of the rocker can be done while the child is standing in the Meerkat.



Read more at R82.co.uk Standing | Meerkat



# Rabbit® Up

The Rabbit Up is a truly adaptable stander. It is secure, reliable and adjustable for all standing needs while also offering the independence to explore using the optional driving wheels. The Rabbit Up enables the child to be an active participant whilst stimulating senses and providing the child with more independence.

## **Highlights about Rabbit Up**



#### **Stander and Mobile stander**

The Rabbit Up not only meets all your standing needs but also offers its users the freedom and opportunity to explore using the two optional driving wheels, giving additional independence to the child.



#### Angle adjustable centre bar

The Rabbit Up has two simple, easy to use angle adjustment options: fixed positions and stepless spindle. This makes transfer easier and helps obtain the best possible position for standing and propelling the driving wheels.



#### **Ouick release wheels**

The two driving wheels are equipped with quick release for easy mounting and dismounting. The wheels are available in a variety of sizes from 22" up to 36" so you can find the perfect fit for the child.



#### **Camber angling**

If you choose to use the driving wheels on your Rabbit Up you have the possibility to change the camber angle. This enables an optimum position of the driving wheels, making propelling easier for the child.

Read more at R82.co.uk Standing | Rabbit Up 121





Mobile stander: Frame, center bar, two driving wheels, standard camber angle, four castors, spoke protectors, fixed positions angle adjustment. Standing frame: Frame, center bar, four castors with brakes.

## **Options, Mobile stander**

#### Four sizes

- Size 1
- Size 2
- Size 3
- Size 4

#### **Driving wheels**

- 22" 24
- 26" 29"
- 32" 36"

## Standard camber angle

- Size 1: 8°
- Size 2+3: 6°
- Size 4: 4°
- Standard plus +3°, +6°, -3°, -6°

#### Push rims

- Open
- Closed

#### Castors

- With brakes
- Without brakes (driving wheels must be fitted)

#### Angle adjustment

- Fixed positions (standard)
- Spindle

## **Options, Standing frame**

#### Four sizes

- Size 1
- Size 2
- Size 3
- Size 4

#### **Brakes**

- 4 castors with brakes (standard)
- Prepared for driving wheels, 4 castors with brakes

#### Angle adjustment

- Fixed positions (standard)
- Spindle

#### Frame colours



Yellow for size 1



Blue for



Anthracite grey for size 3



Silver for



Standard Standing Rabbit Up size 1 (foot plate not standard)



R82 | a part of Etac

Standard Mobile Rabbit Up size 1 (foot plate not standard)



## **Measurements**

		1	2	3	4
<b>Y</b> I	User height: floor to axil	550-740 mm	750-1080 mm	830-1190 mm	960-1420 mm
Ä	Max. load	40 kg	50 kg	60 kg	70 kg
.生.	Angling	0-30°	0-30°	0-30°	0-30°
	Length	700 mm	796 mm	796 mm	926 mm
<u></u>	Length, with driving wheels*	620 mm	750 mm	800 mm	900 mm
<u></u>	Overall width	515 mm	500 mm	600 mm	600 mm
<u>  100  </u>	Overall width, with driving wheels*	635 nn	680 mm	730 mm	720 mm
<b> </b>	Overall height	670-1000 mm	850-1180 mm	990-1310 mm	1200-1530 mm
<u> </u>	Total mass	9 kg	10.5 kg	11 kg	13.5 kg
	Total mass, with driving wheels*	13,5 kg	15,5 kg	16 kg	18,5 kg

Further measurements at our website

123

The Rabbit Up is measured in an upright position, without driving wheels and with the foot plate in lowest position
"With driving wheels: size 1 with 22", size 2 with 26", size 3 with 32", size 4 with 36" wheels

## Nice to know

Rabbit UP offers the child the possibility to independently explore their surroundings by moving around using the optional driving wheels – and maybe even more important to drive away from situations they do not want to participate in.

Read more at R82.co.uk Standing | Rabbit Up





## **Accessories**



Foot plate/Heel stops 86121x/861252



Hinged sandals for flip up foot plate 86122x



Formed foot plates 86149, 86147, 86148



Separate foot support 86120x



Knee supports, hinged 86131x



**Abdominal support** 861411



Hip support, front



Hip support, derotational



Hip support, rear 86150x



Surrounding hip support 861610, 861620, 861630



**Back** support 86401



**Chest support** 86170x

For more accessories - visit our website





Chest support 86180x



Chest support 8618x0



**Distance piece** 861421, 861422



Extension piece 861480



Connecting system for head support 861424-3, 861424-4



Play box 862501



**Tray** 862515

### Gustav, 7 years old

"So you think you can go faster than me..."





125

Read more at R82.co.uk Standing | Rabbit Up



# **Toucan**<sup>®</sup>

The Toucan standing frame comes with soft curves and smart colours. It is designed to support the user in the upright and prone position.

The Toucan is available in 5 sizes.

## **Highlights about Toucan**



Angle adjustable centre bar

The Toucan has easy and smooth angle adjustment making it possible to obtain the optimum standing position.



The Toucan is designed with slots in the centre bar for mounting accessories



#### **Handle with lock**

The handle for angling the centre bar is equipped with a security lock. Pull out the lock whilst turning the handle to change the angle.



#### **Castors with locks**

All four wheels are equipped with locks securing the Toucan for stability. The wheels are very easy to manoeuvre making it easy to move the Toucan around.

Read more at R82.co.uk Standing | Toucan 127





Frame, castors, wheel locks, centre bar, tilt wind system.

# **Options**

#### Five sizes

- Size 1
- Size 2
- Size 3
- Size 4
- Size 5

#### Frame colour







Red for size 1

Blue for size 2+3

Silver for size 4





As standard the Toucan is delivered with frame and centre bar. Please flip once to find accessories

128 Standing | Toucan R82 | a part of Etac



## **Measurements**

		1	2	3	4	5
<b>Y</b> I	User height: floor to axil	340-700 mm	540-900 mm	740-1100 mm	910-1210 mm	1110-1410 mm
	Max. load	30 kg	40 kg	50 kg	60 kg	70 kg
<u></u>	Height floor to frame	60 mm	60 mm	60 mm	60 mm	60 mm
业	Angling	0-30°	0-30°	0-30°	0-20°	0-20°
<u>.</u>	Length	750 mm	880 mm	880 mm	1110 mm	1110 mm
<u>.</u>	Overall width	560 mm	630 mm	630 mm	700 mm	700 mm
	Overall height	670-870 mm	720-1010 mm	970-1230 mm	1040-1340 mm	1240-1540 mm
<u>,                                    </u>	Total mass	7,5 kg	8,5 kg	9 kg	10,5 kg	10,75 kg
					Further measu	rements at our website

## Nice to know

The distance from the floor to the foot plate is only 60 mm, which makes it possible for most of the target users to step into the Toucan.

The target users will be able to maintain antigravity head posture; however, it is possible to mount a head support to compensate for inability to maintain the head alignment during standing.

Read more at R82.co.uk Standing | Toucan





## **Accessories**



Foot plate/Heel stops 86121x/86125x



Hinged sandals 86122x



Formed foot plates 86149, 86147, 86148



**Separate foot support** 86120x



Knee supports 86131x, 86130x



**Abdominal support** 861411



Hip support, front



Hip support, derotational 86161x



Hip support, rear 86150x



**Surrounding Hip support** 8616x0



**Back support** 861401



Chest support 86170x

For more accessories - visit our website







Chest support 86180x



Chest support 8618x0



**Distance piece** 861421, 861422



Extension piece 861480



Connecting system for head **support** 861424-3, 861424-4



86112x



**Tray** 862515



Tray, wood/Arm rests 86110x/ 861111





Read more at R82.co.uk Standing | Toucan







# Caribou<sup>®</sup>

The Caribou is a standing frame that supports a prone or supine standing position. The Caribou is available in 4 sizes and suitable for children and teenagers.

## **Highlights about Caribou**



#### Easy transfer

The Caribou tilts from vertical to horizontal. The horizontal position makes it easier for the carer and less stressful for the user to be placed into position.



#### **Compatible accessories**

The Caribou accessories are compatible with all sizes, and are easily fitted and adjusted on the Caribou frame.



#### **Easy tilt**

The Caribou is supplied with an easy tilt mechanism.



**Anatomic design and comfortable** 

Soft mouldable PU for comfort.

Read more at R82.co.uk Standing | Caribou 133





100 mm castors with locks, height and angle adjustable parted foot support, tilt: From vertical to horizontal.

## **Options**

- Size 1
- Size 2
- Size 3
- Size 4

#### Version\*

- with power - manual

#### Foot plate

Parted foot supports makes it possible to adjust individual foot heights and angles

\*Sizes 3 and 4 are recommended in the power versions to ease tilting

#### Frame colour



Yellow for size 1









Blue for size 3

Anthracite grey for size 4





The Caribou can be ordered as either manual or power version

134 Standing | Caribou R82 | a part of Etac



## **Measurements**

		1	2	3	4
1	User height	850-1120 mm	980-1290 mm	1250-1560 mm	1530-1840 mm
	Max. load	40 kg	60 kg	80 kg	100 kg
<u>,                                    </u>	Height floor to frame	130 mm	130 mm	130 mm	130 mm
业	Angling	0-87°	0-87°	0-87°	0-87°
<u>.</u>	Length	750 mm	920 mm	1110 mm	1110 mm
<u></u>	Overall width	510 mm	570 mm	650 mm	750 mm
	Overall height*	790 mm	910 mm	1130 mm	1410 mm
	Total mass	20,5 kg	24 kg	30,5 kg	35 kg
				Further m	easurements at our websit

\*Configuration: Standard product and top plate (measured from the floor)

## Nice to know

The extensive accessories selection allows children and youngsters who have severe limitations in head and trunk control or who experience limitations in their range of motion to adopt the standing position. Thereby the Caribou supports the child to participate in things they value – and it can bridge the disparities between children with and without disabilities.

When used for prone standing, the shape of the top plate ensures freedom of movement and allows the child to use their arms and hands actively, for instance at the table in front of the body and standing frame.

Read more at R82.co.uk Standing | Caribou





## **Accessories**



Top plate/Head plate 861101-x/861116-x



Cross straps/Chest strap 862102-x/862101-x



Side supports 861102-x



Straps for side supports 86103-x, 862103-x



**Shoulder support** 861126-x



Wedges for side slots 861107-x



**Head support** 9945493-83 / 9967493-83



**Extra cushions** 861106-x



**Back support** 861115-x



**Upper arm supports** 861109-x0, Side supports 861102-x



**Pommel** 861103-x Strap for pommel 862104-x



Knee supports, fixed/ Knee supports, hinged, fixed 861144-x/861144-xx

For more accessories - visit our website





Knee supports, swing-away 861154-x



Knee supports, hinged, swing-away 861154-xx



Formed foot plates 8614x



Heel supports 86125x



**Tray, low mounting** 862131-1, 862131-2



**Tray** 862130-x



Arm rests for plastic tray 861131



Toybar for plastic tray 861135



**Grip handle for plastic tray** 8910750



**Helper's handle** 861108-U



Fix locks 862109-35



Conversion kit 860113-KIT, 860114-KIT




# **Gazelle PS**<sup>®</sup>

Gazelle PS stands for Prone and Supine which means it can be used by children who need to be supported on their fronts or backs. It offers up to 30 degrees of hip abduction in each leg.

The Gazelle PS is available in 3 sizes.

## **Highlights about Gazelle PS**



#### **Tilt vertical-horizontal**

The Gazelle PS tilts from vertical to  $15\,^\circ$  past horizontal. This is an excellent position for respiratory therapy and postural drainage.



#### **Compatible accessories**

The Gazelle accessories are compatible with all sizes, and are easily fitted and adjusted on the Gazelle frame.



#### Multi-adjustable foot supports

Ball joints in foot plates allow for positioning and support for almost any ankle/foot placement.



#### **Manual leg abduction**

Gazelle moves in an anatomically correct motion to help position the head of the femur into the acetabulum. You can adjust either one or both legs to the optimal angle.

Read more at R82.co.uk Standing | Gazelle PS 139





Multi adjustable foot support, knee support, foldable bracket for back support, removable back support, trunk safety belt, castors with lock, manual leg abduction, tilt from vertical to horizontal, adjustable trunk or chest support.

## **Options**

Size - Size 1

- Size 2

- Size 3

#### Frame colour







Yellow for Gi

Green fo

Blue for



As standard the Gazelle PS is supplied with chest and foot straps

## **Measurements**

950 mm 600-105 10 kg 60 k 15° 90:	kg 80 kg
15° 90:	15° 9015°
60 mm 880 m	mm 990 mm
30 mm 580 m	mm 650 mm
60 mm 800 m	mm 900 mm
20 kg 22 k	24 kg
6	60 mm 800 r

140 Standing | Gazelle PS R82 | a part of Etac



## **Accessories**



**Head support system** 88124-xx



Mounting for head support 88130



Facial aperture 88117-x



Anat. head support 88119-x



Side supports 88055-x



Side supports, interior 88075-x



Side supports, exterior 88065-x



Trunk support 88106-x



Angling of lower leg



Strap for side supports 88250-x



Tray with bowl



141

**Tray** 88108-x

For more accessories - visit our website





# Walking, overview

		Mustang	Crocodile	Pony
Target users	Max. load	80 kg	80 kg	50 kg
	User height	700-1300	800-1800	700-1350
Features	Adjustable (forward lean)	•		•
	Height adjustable seat	•		•
	Chest support	•		•
	Handles	•	•	•
	Easy to access	•	•	
	Easy to store and carry	•	•	
	Foldable frame	•	•	
Environment	Indoor	•	•	•
	Outdoor	•*	•	

142 R82 | a part of Etac





# **Mustang**<sup>®</sup>

Mustang is a highly adaptable walking aid for children and youngsters and it can be configured with anterior or posterior support to meet a wide range of walking needs.

The Mustang is available in 4 sizes.

### **Highlights about Mustang**



#### Motivating the child to walk

The user is placed at the centre of the Mustang which provides core stability and enables the user to focus energy on walking.



#### **Easy to access**

The Mustang has been constructed with accessibility in mind. The U-shaped frame and wrap around support, allows for easy access and transfer.



#### **Natural walking pattern**

The Mustang has an easy to adjust central bar for height and angle adjustment. A slight forward leaning position, combined with built-in suspension, stimulates the natural walking pattern.



Accessories are easily added or removed for the carers convenience

Read more at R82.co.uk Walking | Mustang 145





Four sizes, directional stabilizer, drag wheel function, tension adjustable wheels, angle adjustment, reverse brake on/off, height adjustment.

## **Options**

#### Four sizes

- Size 1 yellow
- Size 2 white
- Size 3 red
- Size 4 anthracite grey

#### Wheels

- 8" solid wheels
- 8" slime line wheels
- 4" wheels for size 1 + 2

#### Frame colour









Yellow size 1 White size 2

Red size 3







Standard Mustang size 1

146 Walking | Mustang R82 | a part of Etac



#### **Measurements**

		1	2	3	4
it.	Seat height above floor*	260-440 mm	330-500 mm	450-680 mm	620-840 mm
Ϋ́I	User height: floor to axil*	400-650 mm	600-850 mm	700-1050 mm	870-1300 mm
<u></u>	Max. load	30 kg	40 kg	60 kg	80 kg
<u></u>	Length*	820 mm	900 mm	1000 mm	1100 mm
<u> </u>	Overall width	520 mm	580 mm	650 mm	700 mm
	Turning radius	460 mm	550 mm	600 mm	650 mm
<u>L</u>	Total mass*	13 kg	14 kg	18 kg	19 kg

Further measurements at our website

147

\*Measured with 8" wheels. The Mustang is measured as standard delivered and with the center bar in vertical position. Different settings may result in another set of measurements.

To accomplish some of the measurements above, a slight adjustment might be required. We supply various types of accessories, which is why the weight is listed without

## Nice to know

Mustang is a highly adaptable walking aid, which can be used as a posterior or anterior walker. A recent systematic review was unable to make recommendations for the use of anterior or posterior walkers based on the existing evidence. However, the majority of walking aid user and their parents preferred posterior walkers (Poole 2017).


For children who are "new" walkers the anterior support makes them feel safe and keep them away from hitting objects with the body as they practice their walking skills. The more "experienced" walkers might feel free and more independent and able to get closer to objects with posterior support.

Read more at R82.co.uk Walking | Mustang





## **Accessories**



Chest support 8698x0



Cover 8698x0-x60



Extension piece 869225



Hand bar 86936x



**Abdominal support** 869260



Arm rests, anterior



**Arm prompts** 



Forearm positioner



Sacrum support 86925x



Hip support 8692xx



Hip supports on curved bar 869222 or 869232



**Connecting rod** 869226

For more accessories - visit our website







Sling seat 86957x



**Leg pad** 869470



Guide bar 869510



Storage pouch 869560



Storage bag 869565




Leg separator 8694x0



149

Ankle prompts 8694xx





Read more at R82.co.uk Walking | Mustang



# **Crocodile**®

The Crocodile is a walking aid for children and young adults to help them into a more upright position. It is placed behind the child and its lightweight design makes it easy to move.

The Crocodile is available in 3 sizes.

## **Highlights about Crocodile**



**Height adjustable top frame**Easy to adjust the height without using any tools.



**Multi-adjustable handles**Positioning the grips close to the body gives better support and a better walking position.



#### On/off anti-reverse stops

The anti-reverse wheel stops are used for children with less balance and control.



Foldable frame

Makes storage and transporting the Crocodile easy.

Read more at R82.co.uk Walking | Crocodile 151





Width and height adjustable handles, lightweight durable aluminium frame - sizes 1 and 2, stable steel and aluminium frame - size 3, upholstery and cover for top frame, solid front and rear wheels, swivel front wheels with directional stabilizer, adjustable drag wheels (rear wheels), on/off anti-reverse wheel stop.

## **Options**

#### Three sizes

- Size 1
- Size 2
- Size 3

#### Frame colour



Crocodile size 1

Yellow,

Size 1



Orange, Size 2







Metallic

Anthracite grey, Size 3







Crocodile size 3

152 Walking | Crocodile R82 | a part of Etac



153

### **Measurements**

		1	2	3
4.	User height: floor to handles	420-560 mm	550-690 mm	690-870 mm
	User height, maximum	1150 mm	1500 mm	1800 mm
<u></u>	Max. load	30 kg	45 kg	80 kg
<u>L</u>	Length	720 mm	840 mm	870 mm
, <u> </u>	Overall width	620 mm	690 mm	690 mm
<u>e</u>	Folded length	720 mm	860 mm	1000 mm
	Folded width	620 mm	690 mm	690 mm
i.C.	Folded height	310 mm	270 mm	290 mm
	Turning radius	1000 mm	1140 mm	1140 mm
<u></u>	Total mass	5,5 kg	6,5 kg	12 kg
				Further measurements at our website

## Nice to know

The Crocodile is our hand-held walking aid, allowing children with moderate impairments in standing and walking the possibility to be mobile. The accessories can be used to provide the exact support that is necessary for the user to feel safe and independent

Read more at R82.co.uk Walking | Crocodile





## **Accessories**



Hand brakes 86834/86873



Forearm support 86858/59/60



Forearm positioner 869370



Cross bar 86880



**Back and side support** 



Hip supports 86819



Strap for side supports



**Back support** 86816



Back support, adjustable 86886



Flip down seat 8682x, 86863



Sling seat 8682x



Handle knobs 86861

For more accessories - visit our website











**Bag** 86839

**Bag** Size 3, 86843

**Anti-tips, swing-away** 86815 (1-2), 86815-3S (3)

Attendants handle 86866






155

Read more at R82.co.uk Walking | Crocodile







# **Pony**®

The Pony is a gait trainer for children with disabilities. It has a frame shaped like an arrow and a central bar that holds the chest support and the seat. The Pony is available in 3 sizes.

## **Highlights about Pony**



#### **Back support**

Allows for height and depth adjustment for positioning the pelvis into a neutral position. Furthermore the back support can be removed for easy access into the walker.



**Adjustable forward lean** 

Adjusts to accommodate optimal walking angle and initiate the stepping process.



#### Height adjustable seat

Allows for height adjustment to fit the user and provide optimum positioning when walking. The child is not supposed to sit on the saddle.



#### **Brakes on castors**

The Pony has brakes on the two front castors. Further, it is equipped with one drag wheel and directional stabilizer for added support.

Read more at R82.co.uk Walking | Pony 157





3 sizes, 5 wheels (1 rear and 2 front with lock), bumper wheels, black center plate, hand bar, seat and back support, safety strap.

## **Options**

#### Three sizes

- Size 0 yellow
- Size 1 Purple
- Size 2 Lilac

#### Frame colour








Purple







Standard Pony, purple

Standard Pony, yellow

158 Walking | Pony R82 | a part of Etac



159

## **Measurements**

		0	1	2
rt.	Seat height above floor	240-380 mm	380-500 mm	480-640 mm
<u></u>	Max. load	30 kg	40 kg	50 kg
<u></u>	Length	670 mm	720 mm	1040 mm
<u>.</u>	Overall width	590 mm	690 mm	770 mm
L.	Overall height	520-690 mm	600-820 mm	680-1040 mm
<u></u>	Total mass	7 kg	10 kg	13,5 kg
				Further measurements at our website




Read more at R82.co.uk Walking | Pony





### **Accessories**



Head support system 86424-xx



Helper's handle



Arm rest 8642x



Vertical handles 86432



Sacrum support



Lateral hip supports



Leg protection 86450



Directional stabilizer

For more accessories - visit our website

## Nice to know

The Pony gait trainer is popular in settings, where multiple users, reuse, finically viable and less complex walking solutions are warranted.

The simple and solid design of the Pony assist children, who are unable to maintain a standing position and walk, to be mobile indoor over shorter distances.

160 Walking | Pony R82 | a part of Etac







## **Transportation, overview**

		Stingray Swing-out
Target users	Max load	9/36 kg
	User height	740-1350 mm
Features	Tilt in space	•
	Turnable and swing-out	•
	Lateral support (chest and hips)	•
	Easy-fit system	•
Environment	In car	٠
	Overall width	420 mm



# **Transportation**

Safe and comfortable solutions for the transport of children with special needs!



# Stingray Swing-out® (STSO)

The Stingray Swing-out offers optimum support and function – it tilts, it swivels and it slides. It is available in 4 sizes and suitable for children weighing from 9 to 36 kgs. With its adjustment for growth, it will meet the needs of a growing child for longer.

## **Highlights about Stingray Swing-out**



#### **Swing-out function**

The Swing-out frame enables parents to turn and slide the seat in and out of the car so they don't have to bend when lifting the child into the seat. It also avoids straining their back.



#### **EasyFit System**

The Stingray swing-out car seat frame is easy to fit in cars with ISOFIX but can also be secured using the Easyfit belt tensioning system. The carseat is available in a left or right version.



#### **Tilt Adjustment**

The ingenious tilt has a range of 10-30 degrees making it possible to position the child in an active or more relaxed position.



#### Tested and approved to ECE 4404

The Stingray Swing-out car seat is approved according to ECE 4404 and has been crash tested for user weights between 9 and 36 kgs.





Adjustable tilt in space 10° to 30°, slide and swing-out function on frame, height adjustable 3-point belt guide for the seat belt, height adjustable headrest, adjustable seat depth, easyfit mounting system for Isofix brackets.

## **Options**

4 sizes

- Maat 1

- Maat 2

- Maat 3 - Maat 4

#### Fastening in car

Easyfit system

- Evo tension belt

#### **Cushion colours**

















Frame colour



Black



Standard Stingray Swing-out



## **Measurements**

		1	2	3	4
	Effective seat width	200-280 mm	200-280 mm	260-380 mm	260-380 mm
<u></u>	Effective seat depth	170-240 mm	210-280 mm	270-340 mm	290-360 mm
1	Backrest height	350-450 mm	380-480 mm	430-530 mm	480-580 mm
<u>L</u>	Min. load	9 kg	9 kg	9 kg	9 kg
L	Max load	36 kg	36 kg	36 kg	36 kg
1	Seat plane angle	10-30°	10-30°	10-30°	10-30°
	Overall width	42 cm	42 cm	42 cm	42 cm
	Total mass	21 kg	22 kg	23 kg	24 kg
				Further m	easurements at our websit

## **Accessories**



Side supports NLPE040/041



**Hip supports** 87402x0 (adjustable and non-adjustable)



**5-point belt** ECE Belt NLPE010, 89655-x, 89655-xUK



Outer knee supports NLPE043/044



Foot support 8775010 and NLPE046



Height adjustable headrest



Height adjustable belt guide for proper positioning of the seat belt





## **Toilet and bath, overview**

			Flamingo high-low	Flamingo	Heron
Target users	Max. load		70 kg	70 kg	130 kg
	User height		900-1700 mm	900-1700 mm	1400-1800 mm
Features	Seat plane angle		•	٠	•
	Back recline		•	•	•
	Height adjustable frame		•	•	•
	Foldable frame		•	•	
	Depth adjustable seat				•
Target users	Max. load	Swan	Manatee 75 kg	Orca	Penguin 20 kg
Target users	Max. load User height	Swan  110 kg 1000-1800 mm		Orca 120 kg* 600-1200 mm	
		110 kg	75 kg	120 kg*	20 kg
	User height	110 kg 1000-1800 mm	75 kg	120 kg*	20 kg
	User height Seat plane angle	110 kg 1000-1800 mm	75 kg	120 kg*	20 kg
	User height Seat plane angle Back recline	110 kg 1000-1800 mm	75 kg 700-1300 mm	120 kg*	20 kg
	User height Seat plane angle Back recline Height adjustable frame	110 kg 1000-1800 mm	75 kg 700-1300 mm	120 kg*	20 kg
Target users Features Environment	User height Seat plane angle Back recline Height adjustable frame Foldable frame	110 kg 1000-1800 mm	75 kg 700-1300 mm	120 kg*	20 kg 620-900 mm

168 R82 | a part of Etac



# Toilet & bath

Mobile toilet seats and bath seats provide carers with the best possible working conditions!



# Flamingo High-low®

Flamingo High-low is a modern height adjustable toilet and shower seat, allowing carers to perform bathing tasks in a convenient and good working position. The Flamingo High-low is available in 4 sizes.

## **Highlights about Flamingo High-low**



#### **Adaptation of working conditions**

Height adjustablity eases transition between shower, toilet and sink and provides optimum working conditions for the carer.



#### **Optimum positioning**

The Flamingo High-low offers simple tilt-in-space and back recline functions. It is very easy to position the user with knees higher than the hips which has been proven very effective for bowel function.



#### **Easy to use**

The Flamingo High-low has easily accessible foot pedals for height adjustment and tilt-in-space. The pedals are equipped with security locks for added safety.



Inlay seat with wider sides for comfort and reduced seat width

171

Read more at R82.co.uk Toilet & bath | Flamingo High-low





Height adjustable frame, 4 white castors with lock, tilt in space, foot support, recline, stainless steel frame, active back (no sides), 100 % waterproof upholstery for seat and back, splash guard - high model.

## **Options**

Size

- Size 1

- Size 2

- Size 3 - Size 4

#### **Colours**





Grey

Blue

#### Frame colour



White







Flamingo High-low size 3 and 4



## **Measurements**

Seat		1	2	3	4
	Effective seat width	240 mm	280 mm	320 mm	340 mm
<u></u>	Effective seat depth	250 mm	290 mm	350 mm	400 mm
11_	Backrest height	350 mm	380 mm	420 mm	470 mm
2.	Max. load	35 kg	50 kg	60 kg	70 kg
3	Backrest angle	30°	30°	30°	30°
<u></u>	Length	260 mm	315 mm	360 mm	425 mm
	Overall width	350 mm	380 mm	430 mm	480 mm
11_	Overall height	400 mm	430 mm	480 mm	510 mm
	Total mass	2,5 kg	3 kg	4 kg	4,5 kg
				Further m	neasurements at our websit

Frame and	Seat	1	2	3	4
	Seat surface height at front edge	530-740 mm	530-740 mm	530-740 mm	530-740 mm
<b>1</b>	Max. load	35 kg	50 kg	60 kg	70 kg
之	Seat plane angle	-5°-20°	-5°-20°	-5°-20°	-5°-20°
7	Length	640 mm	640 mm	760 mm	760 mm
	Overall width	500 mm	500 mm	570 mm	570 mm
7	Overall height	860-1080 mm	920-1110 mm	940-1170 mm	1020-1200 mm
7	Total mass	14 kg	14,5 kg	16,5 kg	17 kg

Further measurements at our website





## **Accessories**



**Head support** 990x148-00



Push brace 880110



Sides w/upholstery



Side supports/Tray 83xxx-xxx/88013x



Side Supports 88401-x0



Push handle 880110-HL



Armrests 880152



Commode pan 112401, 880129 Lid for commode pan 8801xx-l



Upholstery with wide sides



Hardware for toilet mounting 88014x



Vest & belts 89xxx-x40



Neoprene hip belt 85423-x

For more accessories - visit our website





Calf Support 880199



**Splash guard, low** 880195, 880196



Potty frame 880160-x



Seat infill 88018x



Suction cups 880150



**Wheels** (100mm) 880100-S



**Back reinforcement kit** 880155



**Back/Rod extension** 880118+880119/880120



Wheels (100 mm) 880100-S





# **Flamingo**<sup>®</sup>

The Flamingo is a toilet/bathing chair. It is adjustable in height and can be placed over the toilet or be used with a pan.

The Flamingo is available in 4 sizes.

## **Highlights about Flamingo**



**Tilt in space and back recline**Tilts and reclines for the child's comfort and positioning.



Accessories can easily be added



**Height adjustable frame**Adjusts to fit the height of the carer.



Adjustable foot supports

Adjust to ensure correct positioning and support.

Read more at R82.co.uk Toilet & bath | Flamingo 177





Height adjustable and foldable frame, 4 white castors with lock, tilt in space, foot support, recline, aluminium frame, active back (no sides), 100 % waterproof upholstery for seat and back, splash guard - high model.

## **Options**

Size

- Size 1

- Size 2 - Size 3

- Size 4

#### **Colours**





Grey

Blue

#### Frame colour



White







R82 | a part of Etac

Standard Flamingo, blue

Toilet & bath | Flamingo



## **Measurements**

Seat		1	2	3	4
	Effective seat width	240 mm	280 mm	320 mm	340 mm
<u> </u>	Effective seat depth	250 mm	290 mm	350 mm	400 mm
11_	Backrest height	350 mm	380 mm	420 mm	470 mm
<b>1</b>	Max. load	35 kg	50 kg	60 kg	70 kg
2	Backrest angle	30°	30°	30°	30°
Ē	Length	260 mm	315 mm	360 mm	425 mm
	Overall width	350 mm	380 mm	430 mm	480 mm
11_	Overall height	400 mm	430 mm	480 mm	510 mm
	Total mass	2,5 kg	3 kg	4 kg	4,5 kg
				Further m	neasurements at our website

Frame an	d Seat	1	2	3	4
	Seat surface height at front edge	500-600 mm	500-600 mm	500-600 mm	500-600 mm
2.	Max. Load	35 kg	50 kg	60 kg	70 kg
之	Seat plane angle	-15°-35°	-15°-35°	-15°-35°	-15°-35°
7	Length	660 mm	660 mm	660 mm	660 mm
	Overall width	440 mm	460 mm	500 mm	540 mm
7.	Overall height	860 mm	890 mm	930 mm	950 mm
0	Folded length	740 mm	740 mm	740 mm	740 mm
	Folded width	440 mm	460 mm	500 mm	540 mm
0	Folded height	200 mm	200 mm	200 mm	200 mm
12.	Total mass	11 kg	13 kg	14 kg	15 kg
				Further m	neasurements at our website

Read more at R82.co.uk

Toilet & bath | Flamingo 179





## **Accessories**



**Head support** 990x148-00



Push brace 880110



Sides w/upholstery



Side supports/Tray 83xxx-xxx/88013x



**Armrests** 880152



Side Supports 88401-x0



Commode pan 112401, 880129 Lid for commode pan



Potty ring 88012x-PR



Hardware for toilet mounting 88014x



Upholstery with wide sides



Vest & belts 89xxx-x40



Neoprene hip belt 85423-x

For more accessories - visit our website





Calf support



Potty frame 880160-x



Seat infill 88018x



**Splash guard, low** 880195, 880196



Suction cups 880150



Wheels (100 mm) 880100-S



Back/Rod extension 880118+880119/880120



**Back reinforcement kit** 880155





Read more at R82.co.uk

Toilet & bath | Flamingo



# **Heron**®

The Heron is a state of the art toilet and bathing seat. It is available with and without power. The simple design and clean appearance conceal multiple functionality features that improve both the user and carer experience.

### **Highlights about Heron**




#### **Optimum positioning**

The Heron is suitable for teenagers and transition to young adults. Change the back height and add a couple of simple accessories to change seat width and seat depth – it is easy to make the Heron fit perfectly!



#### **Perfect working conditions**

One of the great features of the Heron is the hydraulic or powered height adjustment which allows you to adjust the chair into the correct position.



Pressure mapping has been used in the design of Heron to minimise pressure



Heron is available as a hydraulic (100kg) or power version (130kg)

Read more at R82.co.uk

Toilet & bath | Heron 183





### **Standard**

4 grey castors with locks, stainless steel frame, height adjustable frame, integrated push handles, built-in tilt in space, built-in recline, active - height adjustable back, potty ring.

### **Options**

Version

- Hydraulics

- Power

#### **Colours**



Grey

#### Frame colour



White



Heron with hydraulic whithout upholstery (standard)



Heron with power and various accessories

184 Toilet & bath | Heron R82 | a part of Etac



### **Measurements**

		Hydraulics	Power
	Effective seat width	550 mm	550 mm
<u>_</u>	Effective seat depth	450 mm	450 mm
11_	Backrest height	550-620 mm	550-620 mm
	Seat surface height at front edge	510-780 mm	560-780 mm
2.	Max. load	100 kg	130 kg
9	Backrest angle	90°-120°	90°-120°
江	Seat plane angle	-5°-25°	-5°-25°
7	Length	710 mm	710 mm
	Overall width	670 mm	670 mm
7	Overall height	1080-1350 mm	1130-1350 mm
7	Total mass	28 kg	31 kg
			Further measurements at our websi

### Nice to know

Heron is washable at 85°C in a washing machine designed for medical devices. Remember to remove all accessories, unplug all cables and make sure to dismount the control box, battery and remote control before machine washing.

Read more at R82.co.uk

Toilet & bath | Heron







Head support



**Upholstery/Shoulder supports** 880507-04/880508-01



Side supports, swing-away



Arm rests/Tray 880508-03/880508-04



**Hip supports** 880508-09



Commode pan Lid for commode pan




**Upholstery** 101204, back/117005, seat



Splash guard 880508-05



Splash guard, low 880508-07



Seat depth reduction kit 880508-06



Vests & belts 89xxx-x40



Foot supports, short/long 880506-03/880506-01

For more accessories - visit our website










Padding for foot supports 880506-12



Calf support 880509-01



Heel stops 880509-02





Read more at R82.co.uk

Toilet & bath | Heron




# Swan®

Swan is a simple toilet and bathing chair. The seat is made of a smooth material, which makes it easy to clean. The frame and castors are stainless steel with rust-proofed surfaces.

The Swan is available in 6 sizes.

### **Highlights about Swan**



**Tilt in space**Tilts 20° for the child's optimum positioning and comfort.



**Anatomical shape**The Swan has a carefully designed shape that follows the user's contours for improved stability and comfort.



**Holder for pan**For use of pan inside or outside the bathroom.



The stainless steel legs adjust to accommodate most toilets

Read more at R82.co.uk Toilet & bath | Swan 189





### **Standard**

For size 2 - 3: 4 height adjustable legs, holder for commode pan, tilt in space, built-in pommel.

For size 4 - 7: 4 height adjustable legs, holder for commode pan, tilt in space, 4 white castors with lock, removable splash guard.

### **Options**

#### Six sizes

- Size 2
- Size 3
- Size 4
- Size 5
- Size 6
- Size 7

#### Frame and seat colour



White





 $Swan \ is \ available \ in \ 6 \ sizes \ and \ two \ versions. \ Sizes \ 2 - 3 \ include \ a \ built-in \ pommel. \ Sizes \ 4 - 7 \ include \ a \ removable \ splash \ guard.$ 

190 Toilet & bath | Swan R82 | a part of Etac



#### **Measurements**

		2	3	4	5	6	7
	Effective seat width	200 mm	230 mm	250 mm	300 mm	340 mm	410 mm
<u> -</u>	Effective seat depth	270 mm	325 mm	325 mm	380 mm	460 mm	460 mm
11_	Backrest height	380 mm	420 mm	400 mm	480 mm	570 mm	570 mm
2.	Max. load	30 kg	60 kg	70 kg	80 kg	90 kg	110 kg
Ż.	Seat plane angle	18°	18°	27°	22°	20°	20°
7	Length*	460 mm	570 mm	570 mm	570 mm	620 mm	620 mm
<u></u>	Overall width*	490 mm	535 mm	535 mm	535 mm	575 mm	655 mm
7.	Overall height*	680-840 mm	790-950 mm	710-870 mm	780-940 mm	860-1020 mm	860-1020 mm
7	Total mass	8 kg	9,5 kg	8,5 kg	9,5 kg	10,5 kg	11,5 kg

Further measurements at our website 
\*Measured with wheels

191

### **Nice to know**

The simple and easy to clean design makes the swan popular in settings, where multiple users, reuse, finically viable and less complex solutions are warranted. Add the padded upholstery and a head support to increase the comfort of the user.

The swan allows caregivers of children with difficulties in maintaining the sitting position while toileting and showering the possibility to assist them with basic self-care.

Read more at R82.co.uk

Toilet & bath | Swan







**Head support** 990xxxx-00



Push brace 880110



Padded upholstery 8711x



Mesh vest/Cross vest 89xxx-x40



**Tray** 8703x



Commode pan 8801xx Lid for commode pan 8801xx-i



Splash guard 880508-05



**Castors** 87100-20



Rod for foot support 8702x



**Footplate** 87145-x


For more accessories - visit our website



# Cristal, 3 years old "I can't wait to try this chair in the shower"



Read more at R82.co.uk Toilet & bath | Swan



# **Manatee**<sup>®</sup>

The Manatee bathing chair is designed to place the user in exactly the right sitting position. The depth of the seat can be adjusted, allowing the bathing chair to "grow" with the child.

The Manatee is available in 4 sizes.

### **Highlights about Manatee**



**Head cushion/lateral support** 

Provides support and comfort for bathing children. Cushions can be used for head or torso support.



Adjustable back and leg rests

Allows adjustments of the back and lower leg supports with a simple pull of a strap or push of a button.



Provides depth adjustment for growth



195

Adjustable waist strap

Supports the child safely in the Manatee.

Read more at R82.co.uk

Toilet & bath | Manatee





### **Standard**

Head/lateral support, belt for waist/hip.

### **Options**

#### Four sizes

- Size 0
- Size 1
- Size 2
- Size 3

#### **Colours**





Blue

Red

#### Frame colour



Silver



Standard Manatee with Head/lateral support and belt for waist/hip



The Manatee can be lowered to within 14 cm of the bottom of the bath. It is lightweight and it folds flat for storage. The seat can also be placed directly into an ordinary bath tub



### **Measurements**

Seat		0	1	2	3
	Effective seat width	330 mm	330 mm	330 mm	400 mm
<u></u>	Effective seat depth	220-250 mm	300-370 mm	390-460 mm	400-470 mm
11_	Backrest height	430 mm	430 mm	580 mm	640 mm
	Max. load	50 kg	50 kg	60 kg	75 kg
	Overall width	420 mm	420 mm	420 mm	490 mm
=	Folded length	1030-1100 mm	1100-1170 mm	1360-1430 mm	1440-1510 mm
<b>=</b>	Folded width	420 mm	420 mm	420 mm	490 mm
<b>‡</b>	Folded height	130 mm	130 mm	130 mm	130 mm
	Total mass	5,5 kg	5,5 kg	6 kg	7 kg

Frames		Shower frame	Bath tub frame
Ż.	Max. Load	75 kg	75 kg
<u>Z.</u>	Lenght	980 mm	990-1110 mm
<u>П</u> .	Overall width	590 mm	435 mm
7	Overall height*	760-900 mm	270 mm
0	Folded length	980 mm	-
	Folded width	590 mm	-
10-0	Folded height	210 mm	-
<u>.Z.</u>	Total mass	8 kg	2 kg
			*Height measured with chair Further measurements at our website

Read more at R82.co.uk

Toilet & bath | Manatee







**Torso strap** 88302x, 88307x



**Pommel** 88304x, 88309x



**Leg supports** 88305x, 88310x



Pads for head/Lateral support 88301x, 88306x



Bath tub frame - basic version



**Net** 88210x



Shower frame - height adjustable 882100



Wall mounting



For more accessories - visit our website



# Casper, 4 years old 'I look forward to a good long bath'



Read more at R82.co.uk Toilet & bath | Manatee









Max. load 120/20 kg



Environment Indoor only

# Orca | Penguin<sup>®</sup>

The Penguin can be used in the Orca bath tub, but fits ordinary bath tubs as well. With its suction cups, the Penguin sits firmly on most smooth surfaces. Orca and Penguin are available in one size.

### **Highlights about Orca|Penguin**



#### **Excellent working conditions**

Free space around the sides makes it easier for a carer to get close to the child providing optimum comfort. Four smooth running wheels makes it easy to move the bath tub around.



#### **Easy to handle**

The frame can be taken apart to minimise space for storage or transportation. The smooth and hygienic surface makes it very easy to clean.



**Penguin - suction cups** 

Secures the Penguin in the bath tub.



Freedom of movement

Read more at R82.co.uk Toilet & bath | Orca/Penguin 201





# **Standard, Orca**One size, height adjustable frame, soap tray, plug, drain pipe.

# Standard, Penguin One size, suction cups.

#### Frame and seat colour



White





Standard Orca Standard Penguin

202 Toilet & bath | Orca/Penguin R82 | a part of Etac



### **Measurements**

Orca		
	Bath tub width, inside, bottom	400 mm
	Bath tub length, inside, bottom	860 mm
	Max. load	120 kg*
	Length	1110 mm
<u> </u>	Overall width	550 mm
	Overall height	940-1150 mm
	Total mass	16 kg
		*Max user weight: 50 kg

Penguin		
	Inside width, tightest place	160 mm
	Inside length	560-720 mm
	Max. load	20 kg
	Length	800 mm
N	Overall width	330 mm
	Overall height	260 mm
	Total mass	1.35 kg

### **Nice to know**

Orca makes it possible for caregivers using a seated position, such as a wheelchair user, the opportunity to participate and be at eye level in the bathing situation.


The Orca comes with a plug drain and drain pipe for quicker and easier emptying.

Read more at R82.co.uk

Toilet & bath | Orca/Penguin











Extension hose



Orca + Penguin seat



Orca + Flamingo seat



Penguin



Soap dish



Plug













The right accessories are important for obtaining the correct position!







# **Head Supports**

A wide range of head supports provide a large number of opportunities for individual support. Head supports maintain optimum head position and provide support in resting or active positions.

### **Highlights about Head support**



On seating systems

The wide selection of head supports can easily be mounted on seating systems.



On wheelchairs

The wide selection of head supports can easily be mounted on wheelchairs.



**On standing frames** 

The wide selection of head supports can easily be mounted on standing frames.



On bathing chairs

The wide selection of head supports can easily be mounted on bathing chairs.

Read more at R82.com.au Accessories | Head Supports 209







**Head support** no.1 XX=01, no. 2 XX=03, no. 4 xx=06



**Head support** no. 3 xx=05



Head support no. 7, flat XX=13 no. 8, flat, XX=20



Anatomic head support small XX=23, medium. XX=45, large XX=67



Adjustable head support Two laterals XX=77



Adjustable head support with four laterals XX=75



**Anatomic head support** with headband X=5



**Anatomic head support** no. 1, padded. XX=14



Adjustable head support w/ (XX=71) or wo./speakers (XX=73)



Stainless bar system XX=48, XX=47




Swan neck bar system XX=X1, XX=X2




Turtle bar system short xx=92, long xx=93

For more accessories - visit our website









**I-2-I bracket** 87590

Head supports with different shapes and features to suit a wide range of needs



Read more at R82.com.au Accessories | Head Supports



# **Side Supports**

R82 offers a wide range of side supports for seating aids to meet the needs of the individual user. Swing-away side supports are described here – for fixed side supports please see the specific product sections.

### **Highlights about side supports**



For Wombat Living
The wide selection of side supports can easily be mounted on activity chairs.



**For Flamingo**The wide selection of side supports can easily be mounted on bathing chairs.



For Panda Futura and x:panda

The wide selection of side supports can easily be mounted on seating systems.



**For Cougar** 

The wide selection of side supports can easily be mounted on rehabilitation and comfort wheelchairs.

Read more at R82.com.au Accessories | Side Supports 213







**Bars no. 10, 11, 15 and 16**XX=10 (12cm), XX=11 (16cm), XX=15 (Stainless 9,5cm), XX=16
(Stainless 16cm)



**Bar no. 25,26 curved 25mm** xx=25 (length: 12cm) xx=26 (length: 16cm)



**Bar no. 45,46 curved 45 mm** XX=45 (length: 12cm) XX=46 (length: 16cm)



**Bar no. 65,66 curved 65 mm** xx=65 (length: 12cm) xx=66 (length: 16cm)



**Different fittings**Fittings available in multiple shapes for multiple products
83Xxx-xxx



Cushion no. 1, 2 & 3 curved X=1, X=2, X=3



Cushion no. 4 standard



Cushion no. 5 standard



Cushion no. 6 standard



**PU colour**Black, grey and blue 83xxx-xXX
(Blue PU N/A in Australia)



Protection for side supports black = 89150 grey = 89151



**Belt for side support** 84120-25 + 84120-40

For more accessories - visit our website





**Snap lock** 89102-30, 8910702-30

#### **Fitting & Bars**

Different fittings and bar side supports depending on the products.





215

Read more at R82.com.au Accessories | Side Supports



## **Vests and belts**

R82 has developed an extensive range of vests and belts to compliment seating and positioning systems. They increase comfort, safety, and ensure the correct sitting position.

#### **Highlights about Vests and belts**



Support at the hips

A wide selection of hip belts ensure the correct support for all users.



Support at the trunk

A wide selection of trunk supports ensure the correct support for all users.



**Support at the feet** 

Different supports can be added for correct positioning of the user's feet.



**Comfort and protection** 

Read more at R82.com.au Accessories | Vests and belts 217





### **Accessories**



**Shoulder straps** 89663-H



Sternum belt 85430-x



Chest belt 85428-x



**Neoprene vest** 89113-x



Vest 89111-xxx



Cross vest 89625-xxx



**Ankle straps** 85426-x



Belt harness neoprene



Y-hip strap upholstered 85427-x



**Upholstered hip belt** 85424-x



**Neoprene hip belt** 85423-x



**Belt for transportation** 85429-x

For more accessories - visit our website









5-point belt 89665-x



Pelvis harness 89108-x



Hand wrist cuff 85425-x



Foot straps 85xxx



**Neoprene sandals** 85696-1



Thermo bag 87460-xxx



Rain cape with hood 87465-xxx

219





# **Trays**

R82 provide a wide rage of trays with different shapes, colours and materials. Trays with edges to keep toys etc. on the plate, while others are made without to provide more freedom of movement.

### **Highlights about Trays**



Trays with edges

Trays with edges can be mounted on different seating systems.



Swing away trays

Option to choose a tray that can swingaway.



**Transparent tray** 



**Tray with transparent PC** 

Read more at R82.com.au Accessories | Trays 221





## **Accessories**



**Tray with edge, shape 1**Grey ABS 89220-01/ Transparent PET 89230-01



**Tray with edge, shape 2** Wood 89200-2/Grey ABS 89220-2/ PET 89230-2



**Tray with edge, shape 3**Grey ABS 89220-3/ Transparent PET 89230-3



Tray with edge, shape 5 Grey ABS, 8910720-1/Single rod 8910725-1



**Tray with edge, shape 6**Grey ABS 8910720-2, 8910725-2/
8910725-5



Tray with edge, shape 7 Grey ABS 8910725-3



**Tray with edge, shape 8** Available in grey ABS



**Tray without edge, shape 4** Transparent PS 970240-x, 9512x



**Tray without edge, shape 5** Transparent PS 970260-1 and 970260-2



**Tray without edge, shape 6** Transparent PC 87200x0



Toy bar for tray 8910740



**Grip handle for tray** 8910750

For more accessories - visit our website




223



Arm rest for tray 8910730

#### Bertram, 2 years old

"It is easier to reach my toys



Read more at R82.com.au Accessories | Trays



Notes	

Junior   Notes	R82°





#### Find us...

R82 A/S Parallelvej 3 8751 Gedved Denmark Tel +45 79 68 58 88

Fax +45 75 66 51 92 R82@R82.com R82.org

Facebook.com/R82global

Supplier for: Rest of the world

Convaid Products, LLC 2830 California Street Torrance, CA 90503 USA

Tel 844-US-Mobility (844-876-6245)

Convaid-R82.com

facebook.com/ConvaidR82Kids

Etac|R82 Ltd Unit D4A

Coombswood Business Park East Coombswood Way, Halesowen

West Midlands B62 8BH Great Britain

Tel +44 (0)121 561 2222 Fax +44 (0)121 559 5437

R82uk@R82.com R82.co.uk

facebook.com/R82UK

Supplier for:

United Kingdom, Ireland

R82 Inc.

12801 E. Independence Blvd. Matthews, NC 28106

Tel 844-US-Mobility (844-876-6245)

Convaid-R82.com

facebook.com/ConvaidR82Kids

r82@etac.de

R82.de

Etac|R82 GmbH

45699 Herten

Industriestraße 13

Tel. +49(0)2366 500 60

facebook.com/R82.de

Fax +49(0)2366 500 6282

Germany, Luxembourg, Austria,

Supplier for: Switzerland

Supplier for: The Netherlands

R82 au Pty Limited PO Box 284 Botany NSW 1455 Australia

Tel +61 (02) 8213 6666 Fax +61 (02) 8213 6655 R82.com.au

facebook.com/R82Australia au.sales@r82.com

Etac|R82 srl Via Stalingrado 71/7h 40128 Bologna

Etac | R82 Nederland B.V.

Tel +31 (0)513 68 16 86

Fax +31 (0)513 68 19 84

8445 PD Heerenveen

The Netherlands

R82nl@R82.com

Tinweg 8

R82.nl

R82.be

fr.R82.be

Italy

Tel +39 051 2819228 Fax +39 05119873816 Italia.info@r82.com www.R82.it

Supplier for: USA, Canada

Supplier for: USA, Canada Supplier for: Australia

Supplier for: Italy



Read more at R82.co.uk 227



#### Abreast of the future

As the name suggests, R82 was founded in 1982. Since those early pioneering days in the rehabilitation sector, R82 has created simple and user-friendly technical aids and appliances for disabled children and young people.

Flexibility and adaptation to the needs of the individual user were at the heart of the founding of R82 and still is.

Thanks to the consistent philosophy of constant development and high-tech production in a healthy working environment, R82 remains a vibrant, dynamic and growing company.

With a well-designed and high quality range of products in all areas of the technical aid and appliance business, R82 constantly keeps abreast of the future – to the benefit of both users and helpers.



R82 UK Ltd. Unit D4A

Coombswood Business Park
East Coombswood Way
Halesowen
West Midlands B62 8BH
Great Britain
Tel 0121 561 2222
Fax 0121 559 5437
R82uk@R82.com
R82.co.uk